

OPERATING INSTRUCTIONS

Knee walker provides strength and comfort to anyone who may have an injury below the knee. It allows an even distribution of body weight by supporting half of the weight on the device and the other half on the non-injured leg. The walker has a cushion to support both the injured leg and the ankle to make certain that no weight will be placed on the lower leg. Designed for easy steering and stopping. Brakes included.

1. With brake applied, place injured leg on knee pad centered side to side and positioned forward, to cover the full length of the knee pad.
2. With injured leg on knee pad, stand as straight as possible. Adjust the height of the knee walker pad as necessary.
3. With injured foot pointing down, the propulsion leg should be kept as close as possible to the knee pad.
4. Start out with small steps with the propulsion leg, while concentrating on keeping the leg next to the knee pad. With practice, you should be more comfortable with moving forward without deviating to the opposite side.

ASSEMBLY INSTRUCTIONS

1. Remove contents from carton.
2. Insert tiller into front of frame.
3. Raise the tiller and lock in place by sliding the spring loaded lever to the right side (from the knee pad looking forward).
4. Slide locking mechanism into groove on tiller, depress locking handle to secure tiller lock.
5. To adjust the height of the tiller: Slide up or down until it is in a comfortable position. Lock in place with knob and bolt.
6. Insert knee pad post into receptacle on frame. Set to desired height. Secure in place by inserting knob and bolt on left side of frame facing forward and turning clockwise until tight.

HAND BRAKE USE

1. Squeeze the lever to engage the brake.
2. To lock the brake: Pull the lever toward the handle bar and push down the spring loaded push button (located on the top of the brake). When done correctly, push button will remain down and brake lever will be locked in place.
3. To release brake: Squeeze lever. Push button will pop up automatically.

Rhythm Healthcare, LLC
Clearwater, FL 33762

contactus@rhythmhc.com
Website: www.rhythmhc.com
(877) 843-6464

V4: 9/15/23

WARNINGS

- Before using the knee walker, you should be trained by a healthcare professional or dealer.
- Engage the lock and practice bending, reaching and transfers on and off the knee walker.
- The use of the knee walker is individual to each person. You should develop your own methods for use based on your personal level of function and ability.
- DO NOT attempt maneuvers that have not been practiced.
- Be aware of your surroundings. Look for hazards and avoid them.
- Avoid electrical cords, toys, and debris that can be a trip hazard.
- DO NOT use the knee walker near stairs or escalators.
- Use caution when moving from a carpeted area to hard surfaces or when changing levels, such as from a sidewalk to a street.
- DO NOT use knee walker while intoxicated, and/or dizzy.
- Use caution on uneven terrain. DO NOT use on sand.
- Weight cap is 400 lbs.

Limited Lifetime Warranty against manufacturing defects.

This warranty does not cover product failure due to misuse, negligence, improper storage or handling, improper operation or unauthorized modifications.

Normal wear and tear on non-durable components, such as rubber accessories and casters are not covered under this warranty. These items are subject to normal wear and need periodic replacement.

For Warranty Service, it is recommended but not required, that the product be returned to the Rhythm Healthcare, LLC / Lifestyle Mobility Aids dealer through whom it was originally purchased. However, all Rhythm Healthcare, LLC / Lifestyle Mobility Aids dealers are qualified to assist you in obtaining warranty service. If the product is to be returned to Rhythm Healthcare, LLC / Lifestyle Mobility Aids, prior authorization will be required.

In the event of a defect covered by this warranty, Rhythm Healthcare, LLC / Lifestyle Mobility Aids will determine if the product will be repaired or replaced.



ITEM# KN1000

 RHYTHM
HEALTHCARE

KNEE WALKER

400 lb
Weight Capacity

