# P2 PORTABLE OXYGEN CONCENTRATOR



This manual applies to serial numbers 25EBB000001 and later.

## TABLE OF CONTENTS

Symbols.....	1-2
Intended Use.....	3
Contraindications.....	4
General Precautions.....	4-7
Parts Diagram.....	8
User Interface Buttons and Display Symbols.....	9-12
Alerts.....	13-15
Power Supply.....	16-17
Accessories.....	17-18
General Operation.....	19-22
Troubleshooting.....	23
Cleaning the Case.....	24
Cannula Replacement.....	24
Filter Cleaning and Replacement.....	24-25
Battery Care and Maintenance.....	25-26
Disposal of Equipment and Accessories.....	27
Replacement Items.....	27
FAA Patient Travel Regulations.....	27
Specifications.....	28-29
Classifications.....	29
Standards Compliance.....	30
EMC Information.....	31-32
Electromagnetic Emissions.....	32
Warranty.....	33

Symbol	Description
<b>WARNING</b>	A warning indicates that the personal safety of the patient may be involved. Disregarding a warning could result in significant injury
<b>CAUTION</b>	A caution indicates that a precaution or service procedure must be followed. Disregarding a caution could lead to a minor injury or damage to equipment
	See User Manual for Instructions
	AC Power
	DC Power
	U.S. Federal Regulation Restricts this Device to Sale by Order of Physician May also be applicable in other Countries
	No Smoking
	Keep away from open flames
	Keep Dry
	Do not use Oil or Grease
	Do Not Disassemble (contact your equipment provider for servicing by authorized personnel)
	Do Not Dispose of In Unsorted Municipal Waste
	Type BF Applied Part
	Class II (Double Insulated)
	Complies with all applicable EU Directives; including the EU Medical Device Directive
	See Instructions for Use
	Manufacturer
	Authorized representative in the European Community

Symbol	Description
<b>IP22</b>	Protects against solid objects over 12 mm and direct sprays of water up to 15° of vertical
	Date of manufacture
	Serial Number
	This side up
	Fragile
	Storage Humidity (Non-condensing)
	Storage Temperature
	Magnetic Resonance unsafe
	The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft

## **Intended Use**

The P2 Oxygen Concentrator is prescribed for patients requiring supplemental oxygen. It supplies a high concentration of oxygen and is used with a nasal cannula to deliver oxygen from the concentrator to the patient. The P2 is a small and portable device that may be used at the home and can be taken with you while performing your daily activities.

### **WARNING**

This device is not intended to be life-sustaining or life-supporting. This device is not intended for newborn or infant use.

### **WARNING**

A backup oxygen source is recommended for power outages or mechanical problems. Ensure that you have a readily accessible backup oxygen source as advised by your healthcare professional or doctor.

### **CAUTION**

In most countries, this device must be purchased from a doctor or with a doctor's prescription.

### **CAUTION!**

#### **MUST HAVE BACK UP OXYGEN SUPPLY WHEN TRAVELING.**

It is the responsibility of the patient to make back-up arrangements for alternative oxygen supply when traveling. We assume no liability for persons that do not adhere to manufacturer recommendations.

#### **Service Item**

P2 Oxygen Concentrator:  
Molecular Sieve Beds:  
Batteries:

#### **Expected Life**

5 Years  
2 Years  
500 full charge/discharge cycles

### **CAUTION**

The expected life is dependent on the use environment and regular maintenance. Poor conditions will shorten the lifetime of the concentrator.

### **WARNING**

The operator should read and understand this entire manual before using the device.

## **Contraindications**

### **CAUTION**

- This device is not intended to be life sustaining or life supporting.
- Patients who are unable to hear or see an alert from the device, or who are unable to communicate discomfort while wearing the device, will require additional monitoring to avoid injury or harm. If the patient experiences any new symptoms seek medical attention immediately.
- In certain circumstances, oxygen therapy can be hazardous. Please seek medical advice before using this device.
- The P2 is not designed or specified to be used in conjunction with a humidifier, nebulizer or connected with any other equipment. Do not modify the P2 Concentrator. Any modifications performed on the equipment may impair performance or damage equipment and will void your warranty.

## **General Precautions**

### **WARNING**

Oxygen supports combustion. To avoid risk of fire, oxygen therapy should never be used while smoking, while in the same room as someone who is smoking or in the presence of an open flame.

### **WARNING**

Do not submerge the P2 or any of the accessories in liquid.  
Do not expose to water or precipitation.  
Do not operate in the rain/wet weather conditions.  
Exposure to moisture can lead to electrical shock and/or damage.

### **CAUTION**

Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.

### **CAUTION**

Never leave the P2 in high temperatures/high humidity such as in a car in high heat or a bathroom with high humidity. This can damage the device.

### **WARNING**

Geriatric patients or any other patients unable to communicate discomfort, hear or see alarms while using this device, may require additional monitoring.

## **General Precautions-Continued**

### **WARNING**

If you feel any discomfort or are experiencing a medical emergency while using this product, seek medical assistance immediately to avoid harm.

### **WARNING**

Reassess the oxygen delivery settings of this POC periodically to ensure effectiveness of the oxygen therapy.

### **WARNING**

Set the device at the prescribed level and do not increase or decrease your flow rate. Without first consulting with your doctor or healthcare professional.

### **WARNING**

Use this device only as prescribed.  
The use of oxygen therapy can be hazardous in some circumstances.  
Always consult your health care practitioner before using the POC.

### **WARNING**

The settings of the P2 might not correspond with a continuous flow of oxygen.

### **WARNING**

The settings of other models or brands of oxygen therapy equipment do not correspond. With the settings of the P2.

### **WARNING**

There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.

## **General Precautions-Continued**

### **WARNING**

Use only water-based lotions or skin creams that are oxygen compatible during setup or using oxygen therapy. To avoid the risk of fire and burns, never use petroleum or oil-based lotions or salves.

### **WARNING**

Smoking during oxygen therapy is dangerous and may result in fire which can cause serious injury or death of the patient and others.

### **WARNING**

To ensure that you receive the correct therapeutic amount of oxygen delivery according to your medical condition, the P2:

- Must be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
- Must use only the parts and accessories that were provided by the manufacturer and those that were used while your personalized settings were configured.

### **WARNING**

Do not lubricate replaceable fittings, connections, tubing, or other accessories of the oxygen concentrator to avoid the risk of fire and burns.

### **WARNING**

Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.

### **WARNING**

Wind or strong drafts can adversely affect accurate delivery of oxygen therapy.

Examples include:

- Using this equipment beside an open window or in front of a fan.
- Using this equipment in the back seat of an open convertible car.

## **General Precautions-Continued**

### **WARNING**

If any of the following occurs, STOP using immediately and contact your provider:

- Unexplained changes in the performance of this device.
- Unusual or harsh sounds.
- Dropping or mishandling the device or the power supply.
- Water spilled into the enclosure.
- Broken or cracked enclosure.

### **WARNING**

Oxygen is a combustion supporting gas, a fire may start easily if device is used improperly.

Do not leave the nasal cannula on bed coverings or chair cushions if the oxygen concentrator is turned on, but not in use. Always turn the oxygen concentrator off when not in use.

### **WARNING**

To ensure proper function and to avoid the risk of fire and burns:

- Use only with P2 AC and DC power supply.
- Use only with P2 batteries.
- Use only approved P2 accessories.

### **WARNING**

Remove the battery from the device if the P2 will not be used for an extended period.

### **WARNING**

Device operation exceeding the voltage, breath rate, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.

### **WARNING**


Do not modify this system or equipment in any way. Modifications could result in hazards to the user.

**Note: Additional warnings, cautions, and notes are located throughout the manual.**

### **WARNING**

Changes in altitude may affect the amount of oxygen supplied by the device. Consult your physician before travelling to a place with altitude changes.

# Parts Diagram


Flow Setting	Battery Life
1	5 h
2	4 h
3	3 h
4	2 h 14 m
5	1 h 50 m

## User Interface Buttons and Display Symbols


### Display Panel


Symbol	Description	Instructions
	On/Off Button	Press once to turn On. Press and hold for one second to turn Off.
	Audio Alarm Button	<p>Press once to toggle between audible and silent mode.</p> <p>The panel will display the appropriate icon to indicate which mode is enabled:</p> <p>Audible mode----- </p> <p>Silent mode----- </p> <p>When audible mode is enabled, a yellow light will turn on, and a message will display on the screen. Press this button to mute or unmute alarms.</p>
	Flow Setting Control Buttons	Increase or decrease the oxygen flow setting by pressing + or -. Flow settings range from 1 to 5.
	Device Information	Press to display information about the device, including battery temperature, battery status, molecular sieve temperature, molecular sieve runtime, device model, device temperature, device runtime, firmware version, hardware version.

# User Interface Buttons and Display Symbols


## Home Screen

Icon	Description
5	Current flow setting (range is 1 to 5).
	Battery charge level: <ul style="list-style-type: none"><li>• Battery percentage remaining.</li><li>• Battery time remaining (hours: minutes).</li></ul>
	Device runtime since powered on: (hours: minutes).
	Alerts are silenced.
	Audible Alerts are on.

The screen will also display the following icons (see next page):

## User Interface Buttons and Display Symbols

Icon	Description
	Powered by AC or DC only.
	<ol style="list-style-type: none"> <li>1. Powered by battery only (not plugged in and charging).</li> <li>2. Battery level percentage and time remaining (hours: minutes).</li> <li>3. Device On.</li> </ol>
	<ol style="list-style-type: none"> <li>1. Battery is charging.</li> <li>2. Battery level percentage and estimated time to fully charge the battery.</li> <li>3. Device is On.</li> </ol>
	<ol style="list-style-type: none"> <li>1. Battery is charging.</li> <li>2. Battery level percentage and estimated time to fully charge the battery.</li> <li>3. Device is Off.</li> </ol>
	The device has detected an active alarm while in silent mode.
	The device has detected an active alarm while in audible mode.
	Runtime of device (hours: minutes) since it was powered on (Example displays 2 hours 35 minutes).
	Multiple alerts have been detected. Display will scroll to display all alerts.

## User Interface Buttons and Display Symbols

### **Audible Mode:**

Below shows the device has detected an active alarm while in Audible Mode. Example below shows "Absence of Breath".

### **Alerts in Silent Mode:**

The display below shows that the device has detected an active alarm while in Silent Mode. Example below shows "Absence of Breath".

## **Alerts**

### **Adapter plug/unplug:**

An adapter icon displays when adapter is plugged in and disappears when unplugged.  
An audible alarm (if enabled) will also sound.

### **Battery plug/unplug:**

A battery icon displays when battery is connected and disappears when disconnected.  
An audible alarm (if enabled) will also sound.

### **Alarm audio selection:**

An alert will indicate when unit is turned On or Off.

### **Alarm audio pulse duration:**

An audible alert (if enabled) will pulse between 150 ms On, 150 ms off, repeat 2 times.

### **Alarm audio pulse group interval:**

15.5 s (until Alarm returns to normal).

### **Alarm details:**

Reference the table below for additional alarm details.

## Alerts-Continued

<b>Alarm item</b>	<b>Alarm condition</b>	<b>System process</b>	<b>Display of screen</b>
Battery Exhausted	Battery cycle>500 or battery health<50%	Alarm only	Battery Exhausted Replace Battery
Low Input Voltage	Adapter input<17.0 V	Auto-switch to battery until the adapter input>18 V	Low Input Voltage Check Adapter
Absence of Breath	No breath detected continuously for more than 60 seconds	Alarm only	Absence of Breath Check Cannula
Oxygen Concentration <87%	Concentration<87% continuously for more than 10 minutes	Alarm only	Low O2: <87% Contact Provider
Low Battery	5% $\leq$ RSOC $\leq$ 20% Without adapter	Alarm only	Low Battery Charge Now
Oxygen Concentration <50%	Concentration<50% continuously for more than 10 minutes	Alarm only, if improved then alarm clear If condition persists, contact provider	Low O2: <50% Contact Provider
Breath Sensor Fail	Breath sensor failed	Alarm only, if improved then alarm clear If condition persists, contact provider	Breath Sensor Fail Check if Cannula Kinked or Blocked
Oxygen Sensor Fail	Oxygen sensor failed	Alarm only, if improved then alarm clear If condition persists, contact provider	Oxygen Sensor Fail Contact Provider
Gas Delivery Fail	No delivery detected after injection	Alarm only	Gas Delivery Fail Contact Provider
Gas Obstruction	Pipe or nasal blocked	Alarm only	Gas Obstruction Contact Provider
Tank Pressure Fail	Tank pressure failed	Alarm only	Tank Pressure Fail Contact Provider
Sieve Bed Fail	Sieve bed failure or invalid	Auto-shut down after 10 seconds	Sieve Bed Fail Contact Provider

## Alerts-Continued

<b>Alarm item</b>	<b>Alarm condition</b>	<b>System process</b>	<b>Display of screen</b>
Compressor Fail	Compressor failed	Auto-shut down after 10 seconds	Compressor Fail Contact Provider
Valve Check Fail	Valve switch failed	Auto-shut down after 10 seconds	Valve Check Fail Contact Provider
Cooling Fan Fail	Cooling fan failed	Auto-shut down after 10 seconds	Cooling Fan Fail Contact Provider
Battery Depleted	RSOC $\leq$ 5% without adapter	Auto-shut down after 10 seconds	Battery Depleted Replace Battery/ Connect Adapter
System Cold	System temperature $< 32^{\circ}\text{F}/0^{\circ}\text{C}$	Alarm only	System Cold Move to Warmer Place
Battery Cold	Battery temperature $< 32^{\circ}\text{F}/0^{\circ}\text{C}$	Alarm only	Battery Cold Move to Warmer Place
System Hot	System temperature $> 140^{\circ}\text{F}/60^{\circ}\text{C}$	Auto-shut down after 10 seconds	System Hot Move to Cooler Place
Battery Hot	Battery temperature $> 149^{\circ}\text{F}/65^{\circ}\text{C}$	Auto-shut down after 10 seconds	Battery Hot Shut Down, Only Use Adapter
Gas Supply Fail	Flow or concentration below normal after injection	Alarm only	Gas Supply Fail Contact Provider
System Startup Fail	Concentration less than 87% continuously $> 120$ s after system startup	Alarm only	System Startup Fail Contact Provider
Power Supply Fail	System voltage $< 10.5$ V	Auto-shut down after 10 seconds	Power Supply Fail Contact Provider

## Power Supply

### **Standard Lithium-Ion Battery # P2BY001-2**


The P2 is powered by a standard lithium-ion battery. When fully charged, the battery can provide up to 5 hours of operation. Recharge the battery with the AC or DC adapter. Recharging time is not more than 4 hours.

### **AC Power Supply # P2ACA-1**

The AC power supply is used to power the P2 Oxygen Concentrator from an AC power source. When using the AC charger, the power supply automatically adapts to input voltages from 100 V to 240 V (50-60 Hz) allowing it to be used with most power sources throughout the world.

1. Connect A plug to nearest AC power outlet.
2. Connect C to D port.
3. Connect B to the P2.

### **DC Power Supply # P2DCA-1**

The DC power supply is designed for use with the P2 Oxygen Concentrator. The DC power input cable connects directly to the car cigarette lighter or auxiliary DC power supply. The input DC power is 11-16 Vdc, and fuse rating is 15 A/125 V, output is 19 V 6.3 A.

**WARNING** 

Do not use power supplies/adapters or accessories other than those specified above. The use of non-specified accessories may create a safety hazard and/or impair equipment performance.

## **Accessories**

### **Nasal Cannula # P2NC-1**

The P2 Oxygen Concentrator must use a single lumen nasal cannula to provide oxygen to the patient.

**WARNING** 

Nasal cannulas should not be used by more than one person. DO NOT use a cannula that has a length exceeding 25 ft (7.6 m).

**CAUTION**

When using a long cannula, the flow setting may need to be increased. Increasing the cannula length may reduce the perceived noise during oxygen bolus delivery.

**CAUTION**

The nasal cannula is designed for disposable use.

**CAUTION**

Select only FDA approved nasal cannula (e.g. P2NC-1).

### **Carry Bag # P2CB-2**

The P2 carry bag allows you to go out for daily activities.

## **P2 Accessories**

<b>Item Number</b>	<b>Description</b>
P2	Portable Oxygen Concentrator
P2NC-1	Nasal Cannula for the P2/P2-E/P2-E6/P2-E7
P2CB-2	Carry Bag for the P2/P2-E/P2-E6
P2ACA-1	AC Adapter for the P2/P2-E/P2-E6/P2-E7
P2FC-1	Cotton Intake Filter for the P2/P2-E/P2-E6/P2-E7 (pack of 5)
P2BY001-2	Battery for POC (P2/P2-E/P2-E6/P2-E7)
P2DCA-1	DC Adapter for the P2/P2-E/P2-E6/P2-E7

## **General Operation**

### **1. Find a well-ventilated location to place the P2. Ensure that the unit is turned off.**

Ensure the Intake and exhaust has clear access. The P2 should be placed where any alarms can be heard.

#### **WARNING**

Do not use P2 in the presence of flammable anesthetics, detergents, or other chemical vapors.

#### **CAUTION**

Do not block the air intake or air exhaust when operating the equipment. Blocked air circulation or proximity to the heat source can cause internal heat build-up, shut down, or damage to the concentrator.

#### **CAUTION**

The P2 Concentrator is designed for continuous use. It is useful to operate the product frequently for optimal sieve bed life.

#### **CAUTION**

P2 is shipped from factory with battery removed.

**Ensure the Pre-Filter Cover is in place.**

#### **CAUTION**

Do not operate the device without the Cotton Intake Filter and Pre-Filter Cover in place. Operating the device without these in place can damage the device.

## **General Operation-Continued**

### **2. Install the battery.**

Slide the battery into place until the latch returns to the upper position. There will be an audible sound when the battery is in position.

### **3. Connect AC or DC power to P2.**

The green LED on the power adapter will turn on and there will be an audible sound when power is connected.

#### **CAUTION**

Do not place anything in the power supply port other than the supplied wall cord. Avoid the use of electrical extension cords with the P2.

#### **CAUTION**

Power supply is not waterproof. Do not disassemble the power supply.

#### **CAUTION**

When the power is disconnected from the AC outlet, disconnect it from the concentrator to avoid unnecessary battery discharge.

## **General Operation-Continued**

4. Put the P2 into the carry bag, position exhaust vent away from user.

5. Connect the nasal cannula to the nozzle fitting.

Nozzle fitting is located on the top side of the P2 near the Pre-Filter Cover. Connect a nasal cannula to the nozzle fitting on the device. (pictured)

### **CAUTION**

Ensure that the cannula is routed to prevent it from being pinched or kinked to avoid a disruption of oxygen flow.

### **CAUTION**

The cannula is designed for disposable use.

6. Press On/Off  button to turn on P2.

The device will beep, and the indicator light will flash.

"Welcome" will appear on the display at start up. The display will indicate the selected flow setting and power condition. Two-minute warm up time will initiate. During this period the oxygen concentration is building to the specified value but may not have yet reached specification. Under special conditions, a longer warm-up time may be necessary, such as in extremely cold temperatures where the unit was stored or is being operated.

### **CAUTION**

Oxygen concentration may not reach specification during the two-minute warm up time.

### **CAUTION**

30 seconds after startup, the P2 will enter auto-pulse mode. During these 30 seconds, inhalation will not work.

## **General Operation-Continued**

### **7. Set flow rate prescribed by your provider.**

Press the + or - setting buttons to adjust the P2 to the desired flow rate. The current setting can be viewed on the display from 1 to 5.

#### **CAUTION**

Ensure the power is in a well-ventilated place. During operation, the power supply may get hot. Ensure the power supply is cool before handling.

### **8. Wear the nasal cannula on your face and breathe through your nose.**

The P2 will sense if you are breathing from it. If you are not yet breathing through the cannula, the P2 will begin to pulse automatically about every 3 seconds.

As soon as you begin breathing through the cannula, the device will begin delivering pulses based on your breathing. As your breathing rate changes the P2 will sense these changes and adjust the amount of oxygen at your next inhale.

#### **WARNING**

If you feel discomfort using the device, consult your doctor immediately.

#### **CAUTION**

If the oxygen level drops below recommended levels, a Low O2: <87% alert will display on the screen. If the alarm persists, contact your provider.

#### **CAUTION**

If there is no operation for 30 seconds, the display screen will dim. Press any button to activate the display.

#### **CAUTION**

The P2 will notify you with an audible alert (if enabled) and display "Absence of Breath" if no breath has been detected for 60 seconds. After 60 seconds, the device will enter Auto Pulse Mode until breath is detected. Once breath is detected, the device will resume normal oxygen delivery.

#### **General Information**

To disconnect power, unplug the power cord from the wall outlet and disconnect it from the P2.

## Troubleshooting

The table below lists some common problems and solutions. If you can't resolve a problem, please contact your provider.

<b>Problem</b>	<b>Possible Cause</b>	<b>Recommended Solution</b>
Device Won't Turn On	Battery is not installed correctly.	Remove the battery and re-install it correctly.
	Battery is depleted.	Use the AC or DC power adapter operate the device (with the battery inserted) to recharge the battery. If this does not resolve the problem, contact your equipment provider.
	AC supply is not connected properly.	Check AC power connection and verify solid green light on adapter.
	DC Cable is not connected properly.	Check DC power connection at cigarette lighter or auxiliary DC power source.
No Oxygen	The device is not turned on.	Turn on the concentrator.
	Cannula is kinked or obstructed.	Check cannula and its connection to the oxygen outlet port.
	Equipment failure.	Contact your provider.
Oxygen not at full concentration	The device is warming up.	Wait two minutes for the device to warm up. If the issue persists, please contact your provider.
	The sieve beds may require servicing.	Contact your provider to change the sieve beds.
Alarm Occurs	Refer to previous section-Alerts.	Refer to previous section-Alerts.

## **Cleaning the Case**

Clean the outside case with a damp cloth and a solution of mild detergent and water.

### **CAUTION**

Do not allow liquids into any of the controls, the inside or the case, or the oxygen tubing connector. If this occurs, contact your provider for assistance.

### **WARNING**

Do not use alcohol, isopropyl alcohol, ethylene chloride or petroleum-based cleaners on the cases or on the pre-filter covers.

## **Cannula Replacement**

The nasal cannula is disposable. You can buy replacements from your provider and follow the cannula's manufacturer's instructions.

### **CAUTION**

Nasal cannula should be FDA approved and the length should be less than 25 ft (7.6 m).

## **Filter Cleaning and Replacement**

Filters are designed for adequate air flow through the device at the front of the P2.

### **Pre-Filter Cover # P2PFCR-2**

The Pre-Filter Cover must be cleaned weekly to ensure adequate air flow. Clean with a mild liquid detergent and water. Ensure cover is completely dry before use.

### **CAUTION**

It may be necessary to clean the Pre-Filter Cover more often in dusty or polluted environments/conditions.

## Maintenance and cleaning of P2

### **Cotton Intake Filter # P2FC-1**

The Cotton Intake Filter is designed to ensure clean air enters the compressor.

1. Lift Pre-Filter Cover up by bottom end to remove it.
2. Remove Cotton Intake Filter from intake chamber.
3. Install a new Cotton Intake Filter into chamber.
4. Install Pre-Filter Cover.

Pre-Filter Cover and Cotton Intake Filter can be purchased from your provider. In normal conditions, the air filter must be replaced after approximately 3 months of daily use. When subject to conditions with higher levels of dust or dirt, we recommend periodically checking the air filter. If filter is grey or brown color, replace it. The Cotton Intake Filter cannot be cleaned, it must be replaced.

## **Battery Care and Maintenance**

The P2 Lithium-Ion Battery requires special care to ensure proper performance and long life. Only use P2 batteries # P2BY001-2 with your concentrator.

### **CAUTION**

Keep liquids away from the battery. If the battery gets wet, stop using immediately and dispose of battery properly.

### **Battery Replacement**

1. Press down on latch and slide battery out.

2. Insert the P2 battery by sliding battery into place until the latch clicks into place.

## **Battery Care and Maintenance-Continued**

### **Effect of Temperature on Battery Performance**


To extend the run-time of your battery, the device should be used in temperatures between 41°F and 95°F (5°C and 35°C). The number of cycles that the battery will last is highly dependent upon the temperature at which the battery is charged.

### **CAUTION**

We suggest that the room temperature should not exceed 75°F (24°C) when the battery is being charged.

### **Battery Time Remaining Clock**

The P2 continuously displays the battery time remaining. This displayed time is only an estimate and the actual time remaining may vary from this value.

### **CAUTION**

**Store battery in a cool dry place with a charge of 40-50%. BATTERIES SHOULD NOT BE LEFT DORMANT FOR MORE THAN 90 DAYS AT A TIME.**

### **CAUTION**

If the device is not used for an extended period, please remove the battery from the device.

## **Disposal of Equipment and Accessories**

Follow your local governing ordinances for disposal and recycling of the P2 accessories. The battery contains lithium-ion cells and should be recycled and must not be incinerated.

## **P2 Replacement Items**

<b>Part number</b>	<b>Description</b>
P2BY001-2	Battery for POC (P2/P2-E/P2-E6/P2-E7)
P2PFCR-2	Pre-Filter Cover for the P2/P2-E
P2FC-1	Cotton Intake Filter for the P2/P2-E/P2-E6/P2-E7 (pack of 5)

If you need assistance, contact your provider.

## **FAA Patient Travel Regulations**

This device conforms to all applicable FAA acceptance criteria for POC carriage and use onboard aircraft.

### **IMPORTANT**

It is the responsibility of the patient to check with the specific airline carrier when traveling domestically and internationally with a POC. Rhythm Healthcare can assume no liability related to specific airlines rules and regulations.

When traveling with the device, be certain to bring the AC Power Supply and the External Battery Charger (if you have one) with you.

It is advisable to use external power (i.e., plugged into a wall) whenever available to keep the battery fully charged.

Bring enough charged batteries with you to power your concentrator for no less than 150% of the expected duration of your flight, including ground time before and after the flight, security screenings, connections, and a conservative estimate for unanticipated delays.

**Note:** Per FAA regulations, all extra batteries should be individually wrapped and protected to prevent short circuits and carried in carry-on baggage onboard aircraft only.

## System Specifications

### Specifications

Dimensions	L/W/H: 8.70" × 3.35" × 6.30" 22.1 cm × 8.5 cm × 16.0 cm					
Weight	4.37 pounds 1.98 Kg (with battery)					
User Interface	2.8" LCD color display screen					
Sound Level	38 dB(A) (on setting 2)* 1,2					
Time from Turning on the Concentrator to Reach Stated Performance	2 minutes					
Oxygen Concentration	90%-3%/+6% at all settings					
Flow Control Settings and Pulse Volumes	Settings					
		1	2	3	4	5
	Breath Rate	Pulse Volumes(ml)				
	10	21.0	42.0	63.0	84.0	100.0
	15	14.0	28.0	42.0	56.0	66.7
	20	10.5	21.0	31.5	42.0	50.0
	25	8.4	16.8	25.2	33.6	40.0
	30	7.0	14.0	21.0	28.0	33.3
	35	6.0	12.0	18.0	24.0	28.6
	40	5.3	10.5	15.8	21.0	25.0
	±15% at STPD* ±25% over the rated environmental range *STPD is 101.3 kPa at an operating temperature of 68° F, dry					
Breathing Frequency	10 to 40 BPM					
Inspiratory Trigger Sensitivity	-0.12 cm H <sub>2</sub> O					
Time Delay from Onset of Inspiration	<10 ms					
Delivery pressure at the device outlet	Maximum 25 PSI					
Use Mode	Continuous Use					

\* 1 According to ISO 80601-2-69 the Maximum system sound pressure level of the P2 is 42 dB(A) at Setting 2 and 49 dB(A) at Setting 5.

\* 2 Data according to test method 14-1 03/2007 MDS-Hi.

## System Specifications

### Specifications-Continued

Power: AC Power supply DC Power supply Rechargeable Battery	AC Input: 100 to 240 V   AC 50 to 60 Hz DC Input: 11 to 16 V   DC Output: 19 V = 6.3 A Voltage: 14.4 V   DC Rated capacity: 6.8 Ah
Battery Duration	Setting 1: 5 h Setting 2: 4 h Setting 3: 3 h Setting 4: 2 h 14 m Setting 5: 1 h 50 m
Battery Charging Time	Not more than 4 hours.
Environmental Ranges Intended for Operation	Temperature: 41°F to 104°F (5°C to 40°C). Humidity: 10% to 90%, Non-condensing. Altitude: 0 to 10,000 ft (0 to 3048 meters, 70 kPa to 106 kPa).
Environmental Ranges Intended for Shipping And Storage	Temperature: -4°F to 158°F (-20°C to 70°C). Humidity: 5% to 90% Non-condensing. Store in a dry environment. Altitude: 0 to 10,000 ft (0 to 3048 meters, 70 kPa to 106 kPa).
Transportation	Keep Dry. Handle with Care.

### Classifications

Mode of Operation.	Continuous Duty.
Type of Protection Against Electrical Shock.	Class II.
Degree of Protection to Concentrator Components Against Electrical Shock.	Type BF Not intended for cardiac application.
Degree of Protection to Concentrator Components Against Ingress of Water.	IP22-Protects against solid objects over 12 mm and direct sprays of water up to 15° of vertical (IEC 60529).

## **Standards Compliance**

The device is designed to conform to the following standards:

- IEC 60601-1-2, 2014, Medical Electrical Equipment-Part 1-2: General requirements for basic safety and essential performance, Collateral Standard: Electromagnetic disturbances-Requirements and tests.
- IEC 60601-1: Medical Electrical Equipment-Part 1: General Requirements for Basic safety & Essential Performance.
- AAMI ES60601-1: Medical Electrical Equipment-Part 1: General Requirements for Basic safety and Essential Performance.
- IEC 60601-1-8: Medical Electrical Equipment-Part 1-8: General Requirements for Basic Safety and Essential Performance-Collateral Standard: General requirements, tests, and guidance for alarm systems in medical electrical equipment and medical electrical systems.
- IEC 60601-1-11: Medical Electrical Equipment-Part 1-11: General Requirements for Basic Safety and Essential Performance-Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.
- ISO 80601-2-67, Medical Electrical Equipment-Part 2-67: Requirements for basic safety and essential performance of oxygen conserving equipment.
- ISO 80601-2-69, Medical Electrical Equipment-Part 2-69: Requirements for basic safety and essential performance of oxygen concentrator equipment.
- ISO 18562-1: 2017 Biocompatibility, evaluation of breathing gas pathways in healthcare applications-Part 1: Evaluation and testing within a risk management process.
- ISO 18562-2: 2017 Biocompatibility evaluation of breathing gas pathways in healthcare applications-Part 2: Tests for emissions of particulate matter.
- ISO 18562-3: 2017 Biocompatibility evaluation of breathing gas pathways in healthcare applications-Part 3: Tests for emissions of volatile organic compounds (VOCs).
- ISO 10993-1: Biological evaluation of medical devices-Part 1: Evaluation and testing within a risk management process.
- AAMI/ANSIISO 10993-10: 2010, Biological Evaluation of Medical Devices-Part 10: Tests for Skin Irritation.
- AAMI/ANSI/ISO 10993-5: 2009, Biological Evaluation of Medical Devices-Part 5: Tests for in vitro Cytotoxicity.
- ISO 14971: 2019, Medical devices-Application of risk management to medical devices.

## EMC Information

The device has been designed to meet EMC standards throughout its Service Life.

### **Guidance and Manufacturer's Declaration-Electromagnetic Immunity:**

The Concentrator is intended for use in the electromagnetic environment specified below.

The user of the Concentrator should make sure it is used in such an environment.

<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment-Guidance</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV Contact ±15 kV Air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for Power Supply Lines ±1 kV for Input/output Lines	±2 kV for Power Supply Lines ±1 kV for Input output Lines	Main power supply should be that of a typical home or hospital environment.
Surge IEC 61000-4-5	±1 kV Line to Line ±2 V Line to Ground	±1 kV Line to Line ±2 kV Line to Ground	Main power supply should be that of a typical home or hospital environment.
Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	<5% $U_T$ (>95% Dip in $U_T$ ) for 0.5 Cycle at 45 degree increments 70% $U_T$ (30% Dip in $U_T$ ) for 0.5 seconds <5% $U_T$ (>95% Dip in $U_T$ ) for 5 Seconds	<5% $U_T$ (>95% Dip in $U_T$ ) for 0.5 Cycle at 45 degree increments 70% $U_T$ (30% Dip in $U_T$ ) for 0.5 seconds <5% $U_T$ (>95% Dip in $U_T$ ) for 5 Seconds	Main power supply should be that of a typical home or hospital environment. If the user of the Device required continued operation during power outages, it is recommended that the Device be powered from an uninterruptible power supply or battery.
Power Frequency (50/60Hz) Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical home or hospital environment.

Note:  $U_T$  is the A.C. mains voltage prior to application of the test level.

## System Specifications

### **EMC Information-Continued**

Conducted RF	3 Vrms	3 Vrms	<p>Portable and mobile RF communications equipment should not be used closer (to any part of the device, including cables) than the recommended 11.8" 30 cm separation distance.</p> <p>Interference may occur in the vicinity of Equipment marked with the following symbol: </p>
IEC 61000-4-6	150 kHz to 80 MHz	150 kHz to 80 MHz	
Radiated RF	6 Vrms	6 Vrms	
IEC 61000-4-3	Amateur Radio & ISM Bands between 150 kHz and 80 MHz	Amateur Radio & ISM Bands between 150 kHz and 80 MHz	
Radiated RF	10 V/m	10 V/m	
IEC 61000-4-3	80 MHz to 2.7 GHz		

### **Guidance and Manufacturer's Declaration-Electromagnetic Emissions**

The P2 is intended for use in the electromagnetic environment specified below. The user of the P2 should ensure it is used in such an environment.

<b>Emission Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment-Guidance</b>
RF Emissions CISPR 11	Group 1	The P2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The P2 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	

## **Warranty**

The P2 Oxygen Concentrator warranty covers the repair or replacement of the unit. For all units purchased on or after January 1, 2026, the warranty term is 60 months from the date of shipment. For all units purchased before January 1, 2026, the warranty term is 36 months from the date of shipment. Please contact us by telephone or email to return defective equipment under warranty and to resolve any problems. Our trained technicians will help you with any questions or problems with your POC. Please make sure that your returned equipment is packaged safely for transportation, if possible, in its original packaging to avoid damages during shipping. Excluded from the warranty are damages caused by improper usage. Also excluded are replacements of batteries, disposable parts, and consumables. Sieve bed, filters, batteries are expressly excluded from the warranty, except as provided below:

<b>Description</b>	<b>Period</b>
P2 Oxygen Concentrator	5 years - purchased on or after January 1, 2026 3 years - purchased before January 1, 2026
Accessories (Battery, Carry Bag, AC Power Supply and DC Power Supply)	1 year
Sieve Bed	2 years
Disposables (Nasal Cannula, Filters)	No Warranty

Further damage compensation claims of any kind, particularly owing to breach of obligations and unpermitted handling, as well as claims on repayment of expenses paid in vain, are not included in the warranty; the same shall apply to claims on repayments of consequential harm caused by a defect.


Any further claims are excluded in this warranty. The aforementioned limitations do not apply to claims on damages from harm to life, body, or health or attributable to intent or gross negligence, or the product liability law.

This warranty does not cover damage or injury whether to P2 or to personal property or persons caused by accident, misuse, abuse, negligence, failure to install in accordance with Rhythm's installation instructions, failure to operate under conditions of normal use and in accordance with the terms of the operating manual and instructions, failure to maintain in accordance with the applicable service manuals, or alteration or any defects not related to materials or workmanship of P2. This warranty does not cover damage which may occur in shipment. This warranty does not apply to any product or individual part of a product that may have been repaired or altered by anyone other than Rhythm Healthcare or an authorized Rhythm Healthcare service center. This warranty does not apply to any product which is not purchased new.

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Rhythm Healthcare, LLC