TOC3 TRANSPORTABLE OXYGEN CONCENTRATOR


This manual applies to serial numbers 25GBB028195~26?BB??????

TABLE OF CONTENTS

Symbols.....	1-2
Intended Use.....	3
Contraindications.....	4
General Precautions.....	4-8
Parts Diagram.....	9
User Interface Buttons and Display Symbols.....	10-13
Alerts.....	14-15
Built-in Fire Breaker.....	16
Power Supply.....	17-18
Accessories.....	18-21
General Operation.....	22-26
Troubleshooting.....	27
Cleaning the Case.....	28
Cannula Replacement.....	28
Filter Cleaning and Replacement.....	29-31
Battery Care and Maintenance.....	31-32
Disposal of Equipment and Accessories.....	33
Replacement Items.....	33
Specifications.....	34-36
Classifications.....	36
Standards Compliance.....	37
EMC Information.....	38-40
Guidance and Manufacturer's Declaration-Electromagnetic Emissions....	40
Essential Performance List.....	41
Warranty.....	42
Adverse Reactions.....	43

Symbol	Description
WARNING	A warning indicates that the personal safety of the patient may be involved. Disregarding a warning could result in significant injury
CAUTION	A caution indicates that a precaution or service procedure must be followed. Disregarding a caution could lead to a minor injury or damage to equipment
	See User Manual for Instructions
	AC Power
	DC Power
	U.S. Federal Regulation Restricts this Device to Sale by Order of Physician. May also be applicable in other Countries
	No Smoking
	Keep away from open flames
	Keep Dry
	Use No Oil or Grease
	Do Not Disassemble (contact your equipment provider for servicing by authorized personnel)
	Do Not Dispose of In Unsorted Municipal Waste
	Type BF Applied Part
	Class II (Double Insulated)
	See Instructions for Use
	Manufacturer
IP22	Protection against vertically falling water drops
	Date of Manufacture

Symbol	Description
	Serial Number
	This Side Up
	Fragile
	Storage Humidity (Non-condensing)
	Storage Temperature
	Magnetic Resonance unsafe
	The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft

Intended Use

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is intended to provide supplemental low flow oxygen. The device is not intended for life support, nor does it provide patient monitoring capabilities. This device is for adult only.

Intended Population

This device is not intended to be life-sustaining or life-supporting. This device is for adult only.

WARNING

A backup oxygen source is recommended for power outages or mechanical problems. Ensure that you have a readily accessible backup oxygen source as advised by your provider.

CAUTION

In most countries, this device must be purchased from a doctor or with a doctor's prescription.

CAUTION!

MUST HAVE BACK UP OXYGEN SUPPLY WHEN TRAVELING.

It is the responsibility of the patient to make back-up arrangements for alternative oxygen supply when traveling. Rhythm Healthcare assumes no liability for persons that do not adhere to manufacturer's recommendations.

Service Item	Expected Life
Oxygen Concentrator:	5 Years
Sieve Beds:	2 Years
Batteries:	500 full charge/discharge cycles

CAUTION

The expected service life of the concentrator depends on environmental conditions and adherence to regular maintenance. Exposure to poor conditions may significantly reduce the device's lifespan.

WARNING

The operator should read and understand this entire manual before using the device.

Contraindications

- This device is not intended to be life sustaining or life supporting.

WARNING

Your provider has prescribed a specific oxygen flow based on your personal needs. You should not change the flow rate of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator without first consulting with your provider.

CAUTION

- Patients who are unable to hear or see an alert from the device, or are unable to communicate discomfort while using the device, will require additional monitoring to avoid injury or harm. If the patient experiences any new symptoms, seek medical attention immediately.
- In certain circumstances, oxygen therapy can be hazardous. Please seek medical advice before using this device.

General Precautions

WARNING

Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located. If unable to leave the room, you must wait 10 minutes after you have turned the oxygen concentrator off.

WARNING

Do not submerge the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator or any of the accessories in liquid.

Do not expose to water or precipitation.

Do not operate in the rain/wet weather conditions.

Exposure to moisture can lead to electrical shock and/or damage.

CAUTION

Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.

CAUTION

Never leave the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator in high temperatures/high humidity such as in a car in high heat or a bathroom with high humidity. This can damage the device.

General Precautions-Continued

WARNING

If you feel any discomfort or are experiencing a medical emergency while using this product, seek medical assistance immediately to avoid harm.

WARNING

Consult with your provider periodically to reassess the oxygen delivery settings of this the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator to ensure effectiveness of the oxygen therapy.

WARNING

Set the device at the prescribed level and do not increase or decrease your flow rate without first consulting with your provider.

WARNING

Use this device only as prescribed.

The use of oxygen therapy can be hazardous in some circumstances.

Always consult your provider before using the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

WARNING

The settings of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator might not correspond with a continuous flow of oxygen.

WARNING

The settings of other models or brands of oxygen therapy equipment do not correspond with the settings of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

WARNING

There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.

WARNING

Geriatric patients or any other patients unable to communicate discomfort, hear or see alarms while using this device, may require additional monitoring.

General Precautions-Continued

WARNING

Use only water-based lotions or skin creams that are oxygen compatible during setup or using oxygen therapy. To avoid the risk of fire and burns, never use petroleum or oil-based lotions or salves.

WARNING

Smoking during oxygen therapy is dangerous and may result in fire which can cause serious injury or death of the patient and others.

WARNING

To ensure that you receive the correct therapeutic amount of oxygen delivery according to your medical condition, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator:

- Must be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
- Must use only the parts and accessories that were provided by the manufacturer, and those that were used while your personalized settings were configured.

WARNING

Do not lubricate replaceable fittings, connections, tubing, or other accessories of the Oxygen Concentrator to avoid the risk of fire and burns.

WARNING

Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns. Using spare parts not recommended by the manufacturer will void the Warranty.

WARNING

Wind or strong drafts can adversely affect accurate delivery of oxygen therapy.

Examples include:

- Using this equipment beside an open window or in front of a fan.
- Using this equipment in a convertible car with the top down.

WARNING

Use of this device at an altitude above 10,000 ft (3,048 meters) or outside a temperature of 41°F to 104°F (5°C to 40°C) or a relative humidity above 90% is expected to adversely affect the flowrate and the percentage of oxygen and consequently the quality of the therapy.

General Precautions-Continued

WARNING

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will deliver the prescribed flow rate and oxygen concentration reliably after a 5-minute warm-up period.

WARNING

Some respiratory efforts of the patient might not trigger the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

WARNING

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is not intended for use with a tracheotomised patient.

WARNING

Avoid exposure to pollutants and fumes when storing or using the Oxygen Concentrator.

WARNING

The Nasal Cannula is a single patient use medical device. Reusing with other patients may lead to respiratory infections or other illnesses.

WARNING

After powering on the Oxygen Concentrator and allowing it to warm-up for 5 minutes, check for any alarm signals—this includes audible alarms, display alerts, and indicator lights.

- If an alarm is present (excluding the "Absence of Breath" alarm in pulse mode, which clears once breathing is detected), follow the instructions on the main display.
- If the issue persists, discontinue use immediately and contact your provider.
- If no alarms are present, the device is ready for normal operation.

WARNING

To check for air leaks, apply soapy water to the connection between the Nasal Cannula and the oxygen outlet. If bubbles appear, disconnect and reconnect the cannula or contact your provider.

WARNING

If any of the following occurs, STOP using immediately and contact your provider:

- Unexplained changes in the performance of this device.
- Unusual or harsh sounds.
- Dropping or mishandling the device or the power supply.
- Water spilled into the enclosure.
- Broken or cracked enclosure.

General Precautions-Continued

WARNING

Oxygen supports combustion and can cause fires if the device is used improperly.

- Never leave the Nasal Cannula on bed coverings or chair cushions while the Oxygen Concentrator is powered on but not in use.
- Always turn off the Oxygen Concentrator when it is not actively being used.

WARNING

To ensure safe operation and reduce the risk of fire or burns:

- Use only the AC and DC power supplies provided with the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.
- Use only Rhythm Healthcare TOC3 batteries.
- Use only Rhythm Healthcare TOC3 accessories.

WARNING

Remove the battery from the device if the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will not be used for an extended period.

WARNING

Operating the device outside the specified limits for voltage, breath rate, temperature, humidity, or altitude may result in reduced oxygen concentration levels.

WARNING


Do not modify this device or any of its components. Any modifications may create safety hazards and will void the manufacturer's warranty.

WARNING

Altitude changes may impact the oxygen output of the device. Consult your provider before traveling to areas with significant changes in elevation.


Note: Additional warnings, cautions, and notes are located throughout the manual.


Parts Diagram

User Interface Buttons and Display Symbols


Display Panel:


Symbol	Description	Instructions
	On/Off Button	Press once to turn On. Press and hold for three seconds to turn Off.
	Audible Alarm Button	Press once to toggle between Audible and Mute Modes. The panel will display the appropriate icon to indicate which mode is enabled.
	Flow Setting Control Buttons	Adjust the oxygen flow setting by pressing the "+" or "-" buttons.
	Settings Button	Press the button to toggle between Pulse and Continuous Modes. Press and hold to view device information, including battery temperature and status, molecular sieve temperature, runtime, and status, device temperature and runtime, as well as firmware and hardware versions.

User Interface Buttons and Display Symbols-Continued




Home Screen:

Icon	Description
3.0	Current Flowrate (range is 1.0-3.0 in Continuous Mode).
	Battery charge level: <ul style="list-style-type: none"> • Battery percentage remaining. • Battery time remaining (hours:minutes).
	Device runtime since powered on: (hours:minutes).
	Mute Mode
	Audible Mode

The screen will also display the following icons (continued on next page):


User Interface Buttons and Display Symbols-Continued

Icon	Description
	<p>Powered by AC or DC only.</p>
	<ol style="list-style-type: none"> 1. Powered by battery only (not plugged in and charging). 2. Battery level percentage and time remaining (hours: minutes). 3. Device On.
	<ol style="list-style-type: none"> 1. Battery is charging. 2. Battery level percentage and estimated time to fully charge the battery. 3. Device On.
	<ol style="list-style-type: none"> 1. Battery is charging. 2. Battery level percentage and estimated time to fully charge the battery. 3. Device Off.
	<p>The device has detected an active alarm while in Mute Mode.</p>
	<p>The device has detected an active alarm while in Audible Mode.</p>
	<p>Runtime of device (hours:minutes) since it was powered on (Example shows 2 hours 35 minutes).</p>
	<p>Multiple alerts have been detected. Display will scroll to display all alerts.</p>

User Interface Buttons and Display Symbols-Continued


Alerts in Audible Mode:

The device has detected an active alarm in Audible Mode.

Alerts in Mute Mode:

The device has detected an active alarm in Mute Mode.

Alerts

Adapter plug/unplug:

An adapter icon displays when adapter is plugged in and disappears when unplugged. If enabled, an audible alarm will sound.

Battery plug/unplug:

A battery icon displays when battery is connected and disappears when disconnected. If enabled, an audible alarm will sound.

Alarm Item	Alarm Condition	System Process	Display of Screen
Battery Exhausted	Battery cycle > 500 or battery health < 50%	Alarm only	Battery Exhausted Replace Battery
Low Input Voltage	Adapter input < 17.0 V	Auto-switch to battery until the adapter input > 18 V	Low Input Voltage Check Adapter
Absence of Breath	No breath detected continuously for more than 60 seconds	Alarm only	Absence of Breath Check Cannula
Oxygen Concentration < 87%	Concentration < 87% continuously for more than 10 minutes	Alarm only	Low O ₂ : < 87% Contact Provider
Low Battery	5% ≤ RSOC ≤ 20% Without adapter	Alarm only	Low Battery Charge Now
Oxygen Concentration < 50%	Concentration < 50% continuously for more than 10 minutes	If the condition improves, the alarm will clear automatically If the condition persists, contact your provider	Low O ₂ : < 50% Contact Provider
Breath Sensor Fail	Breath sensor failed	If the condition improves, the alarm will clear automatically If the condition persists, contact your provider	Breath Sensor Fail Check if Cannula Kinked or Blocked
Oxygen Sensor Fail	Oxygen sensor failed	If the condition improves, the alarm will clear automatically If the condition persists, contact your provider	Oxygen Sensor Fail Contact Provider

Alerts-Continued

Alarm Item	Alarm Condition	System Process	Display of Screen
Gas Delivery Fail	No delivery detected after injection	Alarm only	Gas Delivery Fail Contact Provider
Gas Obstruction	Pipe or nasal blocked	Alarm only	Gas Obstruction Contact Provider
Tank Pressure Fail	Tank pressure failed	Alarm only	Tank Pressure Fail Contact Provider
Sieve Bed Fail	Sieve bed failure or invalid	Auto-shut down after 10 seconds	Sieve Bed Fail Contact Provider
Compressor Fail	Compressor failed	Auto-shut down after 10 seconds	Compressor Fail Contact Provider
Valve Check Fail	Valve switch failed	Auto-shut down after 10 seconds	Valve Check Fail Contact Provider
Cooling Fan Fail	Cooling fan failed	Auto-shut down after 10 seconds	Cooling Fan Fail Contact Provider
Battery Depleted	RSOC < 5% without adapter	Auto-shut down after 10 seconds	Battery Depleted Replace Battery/ Connect Adapter
System Cold	System temperature < 32°F/0°C	Alarm only	System Cold Move to Warmer Place
Battery Cold	Battery temperature < 32°F/0°C	Alarm only	Battery Cold Move to Warmer Place
System Hot	System temperature > 140°F/60°C	Auto-shut down after 10 seconds	System Hot Move to Cooler Place
Battery Hot	Battery temperature > 149°F/65°C	Auto-shut down after 10 seconds	Battery Hot Shut Down, Only Use Adapter
Gas Supply Fail	Flow or concentration below normal after startup	Alarm only	Gas Supply Fail Contact Provider
Sys Startup Fail	Concentration below 87% continuously > 120 s after system startup	Alarm only	System Startup Fail Contact Provider
Power Supply Fail	System voltage < 10.5 V	Auto-shut down after 10 seconds	Power Supply Fail Contact Provider

Built-in Fire Breaker

Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is with a built-in Fire Breaker that prevents propagation of fire into the device while stopping flow of oxygen. When the device detects that an open flame has spread to the oxygen nozzle through the built-in temperature sensor, the device will stop oxygen flow.

Caution: Only restart the device at least 30 minutes after your detaching of the open flames or fire.

Contact your provider if you have any questions concerning the fire breaker.

Power Supply

Standard Lithium-Ion Battery # TOCBY-2

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is powered by a standard lithium-ion battery.


When fully charged, the TOCBY-2 battery provides up to 6.6 hours of operation. Use the AC or DC Power Supply to recharge the battery; full charging takes no more than 4.5 hours.

AC Power Supply # TOCAC

The AC Power Supply is used to power the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator from an AC power source. It automatically adapts to input voltages ranging from 100 V to 240 V (50-60 Hz), making it compatible with most power sources worldwide.

1. Connect C to D Port.
2. Connect B to the Concentrator.
3. Connect A to the nearest AC power outlet.

DC Power Supply # P2DCA-1

The DC Power Supply is designed for use with the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. The DC power input cable connects directly to a vehicle's cigarette lighter or auxiliary DC power outlet. It accepts an input of 11-16 VDC with a fuse rating of 15 A/125 V, and provides an output of 19 V at 6.3 A.

Power Supply-Continued

WARNING

Use of non-approved the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator accessories may cause damage to the device or user and will void the manufacturer's warranty.

Accessories

Nasal Cannula # TOCNC

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator must requires a single lumen Nasal Cannula to provide oxygen to the patient.

The following Nasal Cannula recommended for use with the Oxygen Concentrator.
Nasal Cannula Model: # TOCNC.
Nasal Cannula Manufacturer: HSINER CO., LTD.

WARNING

Nasal Cannulas should not be used by more than one person. DO NOT use a cannula that has a length exceeding 25 ft (7.6 m).

CAUTION

When using a long cannula, the flow setting may need to be increased. Increasing the cannula length may reduce the perceived noise during oxygen bolus delivery.

CAUTION

The Nasal Cannula is designed for disposable use.

CAUTION

Select FDA approved Nasal Cannula (Rhythm Healthcare Item #: TOCNC).

CAUTION

It is recommended to replace the Nasal Cannula after more than 24 hours of continuous use.

Accessories-Continued

Carry Bag # TOCCB

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Carry Bag can be securely attached to the cart for convenient transport during daily activities. The front pocket is ideal for carrying the AC Power Supply or personal belongings of the patient.

Cart # TOCCRT

Place the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator inside the Carry Bag, then position the Carry Bag on the Cart base. Use the adjustable hook-and-loop fasteners straps to securely fasten the Carry Bag to the Cart handle.

Accessories-Continued

Humidifier (not included)

If your provider has prescribed a humidifier, follow the manufacturer's instructions for use. Attach the humidifier to the oxygen outlet port of the oxygen system. Use of humidifiers not specified for this oxygen system may impair performance and void the warranty.

Select the appropriate humidifier according to the size of the oxygen outlet. The inner diameter of the humidifier connecting pipe should be 4-5 mm.

The following humidifier is recommended for use with the Rhythm Healthcare TOC3 Oxygen Concentrator:

Model: # 30506

Registered Establishment Number: 3003862188

Humidifier Manufacturer: HSINER CO., LTD.

Setting up the Humidifier:

1. Add distilled water to the humidifier. See Figure 1.
2. Secure the humidifier to the device with the strap. See Figure 2.
3. Connect the TOC3 Oxygen Outlet tubing to the Humidifier Oxygen Inlet with the humidifier. See Figure 3.
4. Connect the Nasal Cannula to the Humidifier outlet. See Figure 4.

Figure 1

Figure 2

Figure 3

TOC3 Oxygen
Outlet

Humidifier:
Oxygen Inlet

Figure 4

WARNING

DO NOT use the Humidifier in Pulse Mode.

CAUTION

DO NOT overfill the Humidifier, if used. Fill only with water up to the level indicated by the manufacturer.

CAUTION

The Nasal Cannula must be attached when used with the Humidifier.

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Accessories

Item Number	Description
TOC3	Transportable Oxygen Concentrator
TOCNC	Nasal Cannula
TOCCB	Carry Bag
TOCAC	AC Power Supply
P2DCA-1	DC Power Supply
TOCIF-1	Small Intake Filter
TOCIF-2	Large Intake Filter
TOCBY-2	Battery
TOCCRT	Cart

General Operation

1. Find a well ventilated location to place the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Ensure that the unit is turned off.

Ensure the Intake and exhaust has clear access. The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator should be placed where alarms can be heard.

WARNING

Do not use the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator in the presence of flammable anesthetics, detergents, or other chemical vapors.

CAUTION

The oxygen concentrator requires 30 minutes to warm up from the minimum storage temperature before it is ready for use. Similarly, it requires 30 minutes to cool down from the maximum storage temperature before it is ready for use.

CAUTION

Do not block the Air Intake (Pre-Filter Cover) or Air Exhaust when operating the equipment. Blocked air circulation or proximity to the heat source can cause internal heat build-up, shut down, or damage to the concentrator.

CAUTION

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator can be used 24 hours a day while traveling. Regular use is recommended to preserve optimal performance and extend the life of the sieve beds.

CAUTION

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is shipped from the factory with the battery removed and not installed.

CAUTION


Do not operate the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator without the Small and Large Intake Filters in place. Operating the device without these filters in place can damage the device.


General Operation-Continued

2. Install the battery.

Slide the battery into place until the latch clicks into place. An audible click will confirm that the battery is securely installed.

3. Connect the AC or DC Power Supply to the Rhythm Healthcare TOC3 Transmutable Oxygen Concentrator.

CAUTION

Only use the supplied power cord with the power port. Do not insert any other items into the power port. Do not use extension cords.

CAUTION

Power supply is not waterproof. Do not disassemble the power supply.


CAUTION

When unplugging the power from the wall outlet, also disconnect it from the concentrator to avoid unnecessary battery drain.

General Operation-Continued

4. Connect the Nasal Cannula to the Nozzle Fitting.

Nozzle fitting is located on the top of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Connect the Nasal Cannula to the nozzle fitting on the device (pictured).

CAUTION

Ensure that the cannula is routed to prevent it from being pinched or kinked to avoid a disruption of oxygen flow.

CAUTION

The proper placement and positioning of the prongs of the Nasal Cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.

5. Press the On/Off button to turn on the concentrator.

At startup, the device will beep and the indicator light will flash.

"Welcome" will appear on the display at start up. A five-minute warm-up period will begin, during which the oxygen concentration will gradually increase but may not immediately reach the specified level.

In certain conditions, such as exposure to very cold temperatures during storage or use, a longer warm-up time may be necessary.

CAUTION

Oxygen concentration may not reach specification during the five-minute warm up time.

CAUTION

30 seconds after startup, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will enter Auto-Pulse Mode. Inhalation will not be detected or supported during this initial 30 second period.

CAUTION

Ensure the AC Power Supply is in a well-ventilated area. During operation, the power supply may get hot. Ensure the power supply is cool before handling.

General Operation-Continued

6. Wear the Nasal Cannula on your face and breathe through your nose.


7. Select Pulse or Continuous Mode and adjust flow to the prescribed level.

Press the Settings button to switch between Pulse and Continuous Modes.

Press the "+" or "-" buttons to adjust to your prescribed flow rate.

General Operation-Continued

In Pulse Mode:

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is designed to detect when you are breathing through the cannula. If you haven't started breathing from it yet, the device will automatically deliver a pulse of oxygen approximately every three seconds. Once you begin breathing through the cannula, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will synchronize with your breath, delivering oxygen pulses in response to your inhalation. As your breathing rate changes, the device senses these changes and adjusts the pulse volume accordingly for your next breath.

WARNING

DO NOT use a humidifier in Pulse Mode.

In Continuous Mode:

A constant flow rate of oxygen will be provided.

WARNING

If you feel discomfort using the device, consult your provider immediately.

CAUTION

If the oxygen level drops below recommended levels, a Low O₂: <87% alert will display on the screen. If the alarm persists, contact your provider.

CAUTION

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will notify you with an audible alert (if enabled) and display "Absence Of Breath" if no breath has been detected for 60 seconds. After 60 seconds, the device will enter Auto Pulse Mode until breath is detected. Once breath is detected, the device will resume normal delivery of oxygen.

Note:

To cut off power, unplug the input cord from the AC Power Supply (such as a wall outlet or DC Power Supply), then disconnect it from the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

Troubleshooting

Troubleshooting

The table below lists some common problems and solutions. If you can't resolve a problem, please contact your provider.

Problem	Possible Cause	Recommended Solution
Device Won't Turn On	Battery is not installed correctly.	Remove the battery and re-install it correctly.
	Battery is depleted.	Use the AC or DC Power Supply to operate the device (with the battery inserted) to recharge the battery. If this does not resolve the problem, contact your provider.
	AC Power Supply is not connected properly.	Check AC power connection and verify solid green light on power supply.
	DC Power Supply is not connected properly.	Check DC connection and at cigarette lighter or auxiliary DC power source.
No Oxygen	The device is not turned on.	Turn on the concentrator.
	Cannula is kinked or obstructed.	Check cannula and its connection to the oxygen outlet port.
	Equipment failure.	Contact your provider.
Oxygen Not at Full Concentration	The device is warming up.	Wait 5 minutes for the device to warm-up. If the problem is not solved, please contact your provider.
	The sieve beds may require servicing.	Contact your provider to change the sieve beds.
Alarm Occurs	Refer to previous section-Alerts.	Refer to previous section-Alerts.

Cleaning the Case

The outside case should be cleaned using a damp cloth with a solution of mild detergent and water when necessary.

CAUTION

Do not allow liquids into any of the controls, the interior of the case, or the oxygen tubing connector. If this occurs, contact your provider for assistance.

WARNING

Do not use alcohol, isopropyl alcohol, ethylene chloride or petroleum-based cleaners on the cases or on the filters.

Cannula Replacement

The Nasal Cannula is disposable. Contact your provider should you need Nasal Cannula replacements. Follow the cannula manufacturer's instructions.

CAUTION

Nasal Cannula should be FDA approved and the length should be less than 25 ft (7.6 m).

CAUTION

Nasal Cannula should be rated for 5 liters per minute to ensure proper patient usage and oxygen delivery.

Filter Cleaning and Replacement

Pre-Filter Cover # TOCSP-06

The Pre-Filter Cover is designed for adequate air flow through the device and to ensure clean air is delivered into the compressor.

The Pre-Filter Cover must be cleaned weekly to ensure adequate air flow. Clean with a mild liquid detergent and water. Ensure that the Pre-Filter Cover is completely dry before use.


CAUTION

It may be necessary to clean the Pre-Filter Cover more often in dusty or polluted environments/conditions.

Intake Filters

The filters are designed to allow proper airflow through the bottom of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator includes two intake filters. Clean the filters with a mild liquid detergent and water. Ensure they are completely dry before reinstalling and using them.


Replacing the Small Intake Filter (TOCIF-1)

1. Remove the Filter Box from the bottom of the device.

Filter Cleaning and Replacement-Continued

2. Remove the Small Intake Filter from the Intake Chamber.

3. Insert a new Small Intake Filter into the chamber.
4. Reinstall the Filter Box.

CAUTION

Small Intake Filter (TOCIF-1) should be replaced every 6 months and cleaned once a week. If the device is used in a dusty environment, the Small Intake Filter may need to be replaced more frequently. Do not operate the device without both intake filters installed.

Replacing the Large Intake Filter (TOCIF-2)

1. Remove the Filter Box from the bottom of the device.


2. Remove the Filter Plate.

Filter Cleaning and Replacement-Continued

3. Remove the Large Intake Filter from the Filter Box.

4. Insert a new Large Intake Filter.

5. Place the Filter Plate back in the Filter Box.

6. Reinstall the Filter Box.

CAUTION

The Large Intake Filter (TOCIF-2) should be replaced every year and cleaned once a month. If the device is used in a dusty environment, the filter may need to be replaced more frequently. Do not operate the device without intake filters installed.

Contact your provider for replacement Pre-Filter Covers and Intake Filters.

Battery Care and Maintenance

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Lithium-Ion Battery requires special care to ensure proper performance and long life. Only use the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator # TOCBY-2 battery with your concentrator.

CAUTION

Keep liquids away from battery. If battery gets wet, please stop using immediately and dispose of battery properly.

Battery Replacement

1. Press down on the latch and slide the battery out.
2. Slide battery into place until the latch clicks into place.

Effect of Temperature on Battery Performance

To maximize battery run-time, operate the device in temperatures between 41°F and 95°F (5°C and 35°C). Battery lifespan and performance are strongly influenced by the temperature at which the battery is charged.

CAUTION

The room temperature should not exceed 75°F (24°C) when batteries are being charged.

Battery Time Remaining

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator continuously shows the estimated remaining battery time. Please note that this is only an approximation, and the actual remaining time may differ.

CAUTION

Store the battery in a cool, dry location. To preserve its health and extend its lifespan, charge the battery to 80% every three months if the device will be unused for an extended period. Also, be sure to remove the battery from the device during storage.

BATTERIES SHOULD NOT BE LEFT DORMANT FOR MORE THAN 90 DAYS AT A TIME.

CAUTION

If the device is not used for an extended period of time, please remove the battery from the device and charge it every three months.

Disposal of Equipment and Accessories

Dispose of and recycle Rhythm Healthcare TOC3 Transportable Oxygen Concentrator accessories in accordance with local regulations. The battery contains lithium-ion cells and must be recycled properly-do not incinerate.

Replacement Items

Part number	Description
TOCBY-2	Battery
TOCIF-1	Small Intake Filter
TOCIF-2	Large Intake Filter

If you need assistance, please contact your provider.

System Specifications

Specifications

Mode	Setting											
Continuous Flow	1.0 LPM											
	1.5 LPM											
	2.0 LPM											
	2.5 LPM											
	3.0 LPM											
	Deliverable flow accuracy: The flow error is less than 10% or 0.2 L/min, whichever is greater.											
Mode	Setting	Flowrate L/Min										
Pulse Dose	1	0.2										
	2	0.4										
	3	0.6										
	4	0.8										
	5	1.0										
	6	1.2										
	7	1.4										
	8	1.6										
	9	1.8										
	10	2.0										
Flow Control Settings and Pulse Volumes	Setting	1	2	3	4	5	6	7	8	9	10	
	Breath Rate	Pulse Volumes (mL)										
	10	20	40	60	80	100	120	140	160	180	200	
	15	13	27	40	53	67	80	93	107	120	133	
	20	10	20	30	40	50	60	70	80	90	100	
	25	8	16	24	32	40	48	56	64	72	80	
	30	7	13	20	27	33	40	47	53	60	67	
	35	6	11	17	23	29	34	40	46	51	57	
	40	5	10	15	20	25	30	35	40	45	50	
	Deliverable flow accuracy: $\pm 15\%$ at STPD* $\pm 25\%$ over the rated environmental range *STPD is 101.3 kPa at an operating temperature of 68°F/20°C, dry											
Dimensions (Length x Width x Height)	11.1" L x 6.7" W x 15" H (28.2 cm L x 17.0 cm W x 38.1 cm H)											
Weight	17.2 lbs (7.8 kg) (With TOCBY-2)											

System Specifications

Specifications-Continued

User Interface	2.8" LCD color display screen	
Sound Level	42 dBA at Pulse Setting 2 45 dBA at Continuous Setting 2 LPM (MDS-Hi) ^{*1,2}	
Audio Alarm	≥ 45 dBA	
Time from Turning on the Concentrator to Reach Stated Performance	5 minutes	
Oxygen Concentration	87%-96% (all settings)	
Breathing Frequency	10 to 40 bpm	
Inspiratory Trigger Sensitivity	≤ -4.72" H ₂ O	
Time Delay from Onset of Inspiration	< 10 ms	
Delivery Pressure at The Device Outlet	Pulse Dose: 20 PSI (+10%) Continuous Flow: 5 PSI (+10%)	
Power	AC Power Supply	AC Input: 100 to 240 VAC 50 to 60 Hz
	DC Power Supply	DC Input: 11-16 VDC DC Output: 19 V = 6.3 A
	Rechargeable Battery	TOCBY-2: Voltage: 14.4 VDC Rated capacity: 6.7 Ah + 6.7 Ah
Battery Charging Time	TOCBY-2: Up to 4.5 hours 3.4 Hours at Device Standby Mode when over 90% fully charged by AC power supply. 4.5 Hours at Device Standby Mode when fully charged by AC power supply.	
Environmental Ranges Intended for Operation	Temperature: 41°F to 104°F (5°C to 40°C) Humidity: 10% to 90% Non-condensing Altitude: 0 to 10,000 ft (0 to 3,048 meters, 70 kPa to 106 kPa)	
Environmental Ranges Intended for Shipping and Storage	Temperature: -4°F to 158°F (-20°C to 70°C) Humidity: 5% to 90% Non-condensing Store in a dry environment Altitude: 0 to 10,000 ft (0 to 3,048 meters, 70 kPa to 106 kPa)	

* 1 According to ISO 80601-2-69 the Maximum system sound pressure level of the device is 49 dBA at Pulse Setting 2 and 55 dBA at Continuous Setting 3 LPM.

* 2 Data according to test method 14-1 03/2007 MDS-Hi.

System Specifications

Specifications-Continued

Battery Duration	TOCBY-2	Continuous Flow	1.0 LPM: 3.1 h 1.5 LPM: 2.5 h 2.0 LPM: 1.9 h 2.5 LPM: 1.4 h 3.0 LPM: 1.2 h
		Pulse Dose	Setting 1: 6.6 h Setting 2: 4.9 h Setting 3: 3.9 h Setting 4: 3.6 h Setting 5: 3.3 h Setting 6: 3.0 h Setting 7: 2.6 h Setting 8: 2.3 h Setting 9: 2.0 h Setting 10: 1.8 h
Power Modes	Battery: No Limitation		
	AC Power Supply: No Limitation		
	DC Power Supply: Limited to Settings 1-8 in Pulse Mode Limited to 1-2 LPM in Continuous Mode		
Operating Time	24-hours a day when connected to an external AC or DC power source.		
Transportation	Keep Dry, Handle With Care.		

Classifications

Mode of Operation	Continuous Duty
Type of Protection Against Electrical Shock	Class II
Degree of Protection to Concentrator Components Against Electrical Shock	Type BF Not intended for cardiac application
Degree of Protection to Concentrator Components Against Ingress of Water	IP22-Protects against solid objects over 12.5 mm and direct sprays of water up to 15° of vertical (IEC 60529)

Standards Compliance

The device is designed to conform to the following standards:

- IEC 60601-1: 2005/ADM1: 2012/ADM2: 2020, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2: 2014/AMD1: 2020, Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances-Requirements and tests
- IEC 60601-1-6: 2010/AMD1: 2013/AMD2: 2020, Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
- IEC 60601-1-8: 2006/AMD1: 2012/AMD2: 2020, Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
- IEC 60601-1-11: 2015/AMD1: 2020, Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- ISO 80601-2-67: 2020, Medical electrical equipment - Part 2-67: Particular requirements for basic safety and essential performance of oxygen - conserving equipment
- ISO 80601-2-69: 2020, Medical electrical equipment - Part 2-69: Particular requirements for the basic safety and essential performance of oxygen - concentrator equipment
- ISO 18562-1: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 1: Evaluation and testing within a risk management process
- ISO 18562-2: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter
- ISO 18562-3: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs)
- ISO 10993-1: 2018, Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- ISO 10993-5: 2009, Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- ISO 10993-10: 2021, Biological evaluation of medical devices - Part 10: Tests for skin sensitization
- ISO 10993-23: 2021, Biological evaluation of medical devices - Part 23: Tests for irritation

EMC Information

The device has been designed to meet EMC standards throughout its Service Life.

Guidance and Manufacturer's Declaration-Electromagnetic Immunity:

The Concentrator is intended for use in the electromagnetic environment specified below.

The users of the Concentrator should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air	± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	± 2 kV for Power Supply Lines ± 1 kV for Input/ Output Lines	± 2 kV for Power Supply Lines ± 1 kV for Input/ Output Lines	Mains power quality should be that of a typical home or hospital environment.
Surge IEC 61000-4-5	± 1 kV Line to Line	± 1 kV Line to Line	Main power quality should be that of a typical home or hospital environment.
Voltage Dips, Short Interruptions on Power Supply Input Lines IEC 61000-4-11	< 5% U_T (> 95% Dip in U_T) for 0.5 Cycle at 45 degree increments 70% U_T (30% Dip in U_T) for 0.5 Seconds < 5% U_T (> 95% Dip in U_T) for 5 Seconds	< 5% U_T (> 95% Dip in U_T) for 0.5 Cycle at 45 degree increments 70% U_T (30% Dip in U_T) for 0.5 Seconds < 5% U_T (> 95% Dip in U_T) for 5 Seconds	Mains power quality should be that of a typical home or hospital environment. If the users of the device be required continued operation during power mains interruptions, it is recommended that device be powered from an uninterruptible power supply or battery.
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical home or hospital environment.
Note: U_T is the A.C. mains voltage prior to application of the test level.			

EMC Information-Continued

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended 12 inches separation distance. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	6 Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz	6 Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7G Hz	10 V/m	

CAUTION

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is may be used in home, institution and various mobile environments.

CAUTION

Electromagnetic interference could disturb performance of the device (oxygen concentration and flow), resulting in an alarm.

When used in the above environment, it is also necessary to avoid exposure to known sources of Electro-Magnetic Interference, including MRI, RFID, diathermy, and electrocautery.

CAUTION

Use of accessories (Battery, AC Power Supply, DC Power Supply) other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic, immunity of this equipment and result in improper operation.

CAUTION

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches to any part of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator, including cables supplied by the manufacturer. The use of this device closer that the recommended distance can result in a degradation of the performance of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

EMC Information-Continued

CAUTION

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Guidance and Manufacturer's Declaration-Electromagnetic Emissions

The Concentrator is intended for use in the electromagnetic environment specified below. The user of the Concentrator should make sure it is used in such an environment.

Emissions test	Compliance	Electromagnetic Environment-Guidance
RF Emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The device is suitable for uses in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	

Essential Performance List

Essential Performances	Relevant Requirements or Conditions
Accuracy of concentration	Within 87% - 96% (at all settings)
Accuracy of oxygen delivery	Within 25% as indicated in the instruction for use
Power Supply Failure Alert	System voltage < 10.5 V
Internal electrical power source nears depletion	RSOC < 5% without adapter
Low Oxygen Concentration Alert	The oxygen concentration remains below 87% continuously for more than 600 seconds
Start Up Conditions Alert	The oxygen concentration remains below 87% continuously for more than 120 seconds after system startup
Triggering signal absence alarm condition (Only in Pulse Mode)	No breath has been detected continuously for more than 60 seconds.
Gas Supply Failure Alert	Flow or concentration below normal after injection

Warranty

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator warranty covers the repair or replacement of the unit. The warranty term is 36 months from the date of shipment. Please contact us by telephone or email to return defective equipment under warranty and to resolve any problems. Our trained technicians will help you with any questions or problems with your the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Please make sure that your returned equipment is packaged safely for transportation, if possible, in its original packaging to avoid damages during shipping. Excluded from the warranty are damages caused by improper usage. Also excluded are replacements of batteries, disposable parts, and consumables. Sieve bed, filters, batteries are expressly excluded from the 36 months warranty, except as provided below:

Description	Period
Oxygen Concentrator	3 Years
Accessories (Battery, Carry Bag, AC Power Supply and DC Power Supply)	1 Year
Disposables (Nasal Cannula, Filters)	No Warranty
Sieve Bed	2 years - purchased on or after January 1, 2026 3 years - purchased before January 1, 2026

Further damage compensation claims of any kind, particularly owing to breach of obligations and unpermitted handling, as well as claims on repayment of expenses paid in vain, are not included in the warranty; the same shall apply to claims on repayments of consequential harm caused by a defect.

Any further claims are excluded in this warranty. The aforementioned limitations do not apply to claims on damages from harm to life, body, or health or attributable to intent or gross negligence, or the product liability law.

This warranty does not cover damage or injury whether to the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator or to personal property or persons caused by accident, misuse, abuse, negligence, failure to install in accordance with Rhythm's installation instructions, failure to operate under conditions of normal use and in accordance with the terms of the operating manual and instructions, failure to maintain in accordance with the applicable service manuals, or alteration or any defects not related to materials or workmanship of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. This warranty does not cover damage which may occur in shipment. This warranty does not apply to any product or individual part of a product that may have been repaired or altered by anyone other than Rhythm Healthcare or an authorized Rhythm Healthcare service center. This warranty does not apply to any product which is not purchased new.

Adverse Reactions

MedWatch is the Food and Drug Administration's (FDA) program for reporting serious reactions, product quality problems, therapeutic inequivalence/failure, and product use errors with human medical products, including drugs, biologic products, medical devices, dietary supplements, infant formula, and cosmetics.

If you think you or someone in your family has experienced a serious reaction to a medical product, you are encouraged to take the reporting form to your doctor. Your health care provider can provide clinical information based on your medical record that can help FDA evaluate your report.

However, we understand that for a variety of reasons, you may not wish to have the form filled out by your health care provider, or your health care provider may choose not to complete the form. Your health care provider is not required to report to the FDA. In these situations, you may complete the Online Reporting Form yourself.

You will receive an acknowledgement from FDA when your report is received. Reports are reviewed by FDA staff. You will be personally contacted only if we need additional information.

Submitting Adverse Event Reports to FDA

Use one of the methods below to submit voluntary adverse event reports to the FDA at <http://www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home>

Consumer Reporting Form FDA 3500B.

Follow the instructions on the form to either fax or mail it in for submission. For help filling out the form, see MedWatchLearn.

The form is available at <http://www.fda.gov/downloads/aboutFDA/reportsmanuals/forms/forms/ucm349464.pdf>.


Call FDA at 1-800-FDA-1088 to report by telephone. Reporting Form FDA 3500 is commonly used by health professionals.

The form is available at <http://www.fda.gov/downloads/aboutFDA/reportmanuals/forms/forms/ucm163919.pdf>.

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Rhythm Healthcare, LLC