# TOC3 TRANSPORTABLE OXYGEN CONCENTRATOR

# TABLE OF CONTENTS

Symbols.....1-2

Intended Use.....3

Contraindications.....4

General Precautions.....4-8

Parts Diagram.....9

User Interface Buttons and Display Symbols.....10-13

Alerts.....14-15

Power Supply.....16-17

Accessories.....17-20

General Operation.....21-25

Troubleshooting.....26

Cleaning the Case.....27

Cannula Replacement.....27

Filter Cleaning and Replacement.....28-30

Battery Care and Maintenance.....30-31

Disposal of Equipment and Accessories.....32

Replacement Items.....32

Specifications.....33-35

Classifications.....35

Standards Compliance.....36

EMC Information.....37-39


Guidance and Manufacturer's Declaration-Electromagnetic Emissions....39

Essential Performance List.....40

Warranty.....41

Adverse Reactions.....42

| Symbol                                                                              | Description                                                                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>WARNING</b>                                                                      | A warning indicates that the personal safety of the patient may be involved. Disregarding a warning could result in significant injury                  |
| <b>CAUTION</b>                                                                      | A caution indicates that a precaution or service procedure must be followed. Disregarding a caution could lead to a minor injury or damage to equipment |
|    | See User Manual for Instructions                                                                                                                        |
|    | AC Power                                                                                                                                                |
|    | DC Power                                                                                                                                                |
|    | U.S. Federal Regulation Restricts this Device to Sale by Order of Physician. May also be applicable in other Countries                                  |
|    | No Smoking                                                                                                                                              |
|    | Keep away from open flames                                                                                                                              |
|    | Keep Dry                                                                                                                                                |
|    | Use No Oil or Grease                                                                                                                                    |
|    | Do Not Disassemble (contact your equipment provider for servicing by authorized personnel)                                                              |
|    | Do Not Dispose of In Unsorted Municipal Waste                                                                                                           |
|  | Type BF Applied Part                                                                                                                                    |
|  | Class II (Double Insulated)                                                                                                                             |
|  | See Instructions for Use                                                                                                                                |
|  | Manufacturer                                                                                                                                            |
| <b>IP22</b>                                                                         | Protection against vertically falling water drops                                                                                                       |
|  | Date of Manufacture                                                                                                                                     |

| Symbol                                                                            | Description                                                                                                                                    |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Serial Number                                                                                                                                  |
|  | This Side Up                                                                                                                                   |
|  | Fragile                                                                                                                                        |
|  | Storage Humidity (Non-condensing)                                                                                                              |
|  | Storage Temperature                                                                                                                            |
|  | Magnetic Resonance unsafe                                                                                                                      |
|  | The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft |

## **Intended Use**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is intended to provide supplemental low flow oxygen. The device is not intended for life support, nor does it provide patient monitoring capabilities. This device is for adult only.

### **Intended Population**

This device is not intended to be life-sustaining or life-supporting. This device is for adult only.

### **WARNING**

A backup oxygen source is recommended for power outages or mechanical problems. Ensure that you have a readily accessible backup oxygen source as advised by your provider.

### **CAUTION**

In most countries, this device must be purchased from a doctor or with a doctor's prescription.

### **CAUTION!**

#### **MUST HAVE BACK UP OXYGEN SUPPLY WHEN TRAVELING.**

It is the responsibility of the patient to make back-up arrangements for alternative oxygen supply when traveling. Rhythm Healthcare assumes no liability for persons that do not adhere to manufacturer's recommendations.

| <b>Service Item</b>  | <b>Expected Life</b>             |
|----------------------|----------------------------------|
| Oxygen Concentrator: | 5 Years                          |
| Sieve Beds:          | 2 Years                          |
| Batteries:           | 500 full charge/discharge cycles |

### **CAUTION**

The expected service life of the concentrator depends on environmental conditions and adherence to regular maintenance. Exposure to poor conditions may significantly reduce the device's lifespan.

### **WARNING**

The operator should read and understand this entire manual before using the device.

## **Contraindications**

- This device is not intended to be life sustaining or life supporting.

### **WARNING**

Your provider has prescribed a specific oxygen flow based on your personal needs. You should not change the flow rate of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator without first consulting with your provider.

### **CAUTION**

- Patients who are unable to hear or see an alert from the device, or are unable to communicate discomfort while using the device, will require additional monitoring to avoid injury or harm. If the patient experiences any new symptoms, seek medical attention immediately.
- In certain circumstances, oxygen therapy can be hazardous. Please seek medical advice before using this device.

## **General Precautions**

### **WARNING**

Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located. If unable to leave the room, you must wait 10 minutes after you have turned the oxygen concentrator off.

### **WARNING**

Do not submerge the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator or any of the accessories in liquid.

Do not expose to water or precipitation.

Do not operate in the rain/wet weather conditions.

Exposure to moisture can lead to electrical shock and/or damage.

### **CAUTION**

Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.

### **CAUTION**

Never leave the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator in high temperatures/high humidity such as in a car in high heat or a bathroom with high humidity. This can damage the device.

## **General Precautions-Continued**

### **WARNING**

If you feel any discomfort or are experiencing a medical emergency while using this product, seek medical assistance immediately to avoid harm.

### **WARNING**

Consult with your provider periodically to reassess the oxygen delivery settings of this the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator to ensure effectiveness of the oxygen therapy.

### **WARNING**

Set the device at the prescribed level and do not increase or decrease your flow rate without first consulting with your provider.

### **WARNING**

Use this device only as prescribed.  
The use of oxygen therapy can be hazardous in some circumstances.  
Always consult your provider before using the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

### **WARNING**

The settings of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator might not correspond with a continuous flow of oxygen.

### **WARNING**

The settings of other models or brands of oxygen therapy equipment do not correspond with the settings of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

### **WARNING**

There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.

### **WARNING**

Geriatric patients or any other patients unable to communicate discomfort, hear or see alarms while using this device, may require additional monitoring.

## **General Precautions-Continued**

### **WARNING**

Use only water-based lotions or skin creams that are oxygen compatible during setup or using oxygen therapy. To avoid the risk of fire and burns, never use petroleum or oil-based lotions or salves.

### **WARNING**

Smoking during oxygen therapy is dangerous and may result in fire which can cause serious injury or death of the patient and others.

### **WARNING**

To ensure that you receive the correct therapeutic amount of oxygen delivery according to your medical condition, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator:

- Must be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
- Must use only the parts and accessories that were provided by the manufacturer, and those that were used while your personalized settings were configured.

### **WARNING**

Do not lubricate replaceable fittings, connections, tubing, or other accessories of the Oxygen Concentrator to avoid the risk of fire and burns.

### **WARNING**

Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns. Using spare parts not recommended by the manufacturer will void the Warranty.

### **WARNING**

Wind or strong drafts can adversely affect accurate delivery of oxygen therapy.

Examples include:

- Using this equipment beside an open window or in front of a fan.
- Using this equipment in a convertible car with the top down.

### **WARNING**

Use of this device at an altitude above 10,000 ft (3,048 meters) or outside a temperature of 41°F to 104°F (5°C to 40°C) or a relative humidity above 90% is expected to adversely affect the flowrate and the percentage of oxygen and consequently the quality of the therapy.

## **General Precautions-Continued**

### **WARNING**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will deliver the prescribed flow rate and oxygen concentration reliably after a 5-minute warm-up period.

### **WARNING**

Some respiratory efforts of the patient might not trigger the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

### **WARNING**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is not intended for use with a tracheotomised patient.

### **WARNING**

Avoid exposure to pollutants and fumes when storing or using the Oxygen Concentrator.

### **WARNING**

The Nasal Cannula is a single patient use medical device. Reusing with other patients may lead to respiratory infections or other illnesses.

### **WARNING**

After powering on the Oxygen Concentrator and allowing it to warm-up for 5 minutes, check for any alarm signals—this includes audible alarms, display alerts, and indicator lights.

- If an alarm is present (excluding the "Absence of Breath" alarm in pulse mode, which clears once breathing is detected), follow the instructions on the main display.
- If the issue persists, discontinue use immediately and contact your provider.
- If no alarms are present, the device is ready for normal operation.

### **WARNING**

To check for air leaks, apply soapy water to the connection between the Nasal Cannula and the oxygen outlet. If bubbles appear, disconnect and reconnect the cannula or contact your provider.

### **WARNING**

If any of the following occurs, STOP using immediately and contact your provider:

- Unexplained changes in the performance of this device.
- Unusual or harsh sounds.
- Dropping or mishandling the device or the power supply.
- Water spilled into the enclosure.
- Broken or cracked enclosure.

## **General Precautions-Continued**

### **WARNING**

Oxygen supports combustion and can cause fires if the device is used improperly.

- Never leave the Nasal Cannula on bed coverings or chair cushions while the Oxygen Concentrator is powered on but not in use.
- Always turn off the Oxygen Concentrator when it is not actively being used.

### **WARNING**

To ensure safe operation and reduce the risk of fire or burns:

- Use only the AC and DC power supplies provided with the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.
- Use only Rhythm Healthcare TOC3 batteries.
- Use only Rhythm Healthcare TOC3 accessories.

### **WARNING**

Remove the battery from the device if the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will not be used for an extended period.

### **WARNING**

Operating the device outside the specified limits for voltage, breath rate, temperature, humidity, or altitude may result in reduced oxygen concentration levels.

### **WARNING**


Do not modify this device or any of its components. Any modifications may create safety hazards and will void the manufacturer's warranty.

### **WARNING**

Altitude changes may impact the oxygen output of the device. Consult your provider before traveling to areas with significant changes in elevation.


**Note: Additional warnings, cautions, and notes are located throughout the manual.**

## Parts Diagram

## User Interface Buttons and Display Symbols


### Display Panel:


| Symbol | Description                  | Instructions                                                                                                                                                                                                                                                                        |
|--------|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|        | On/Off Button                | Press once to turn On.<br>Press and hold for three seconds to turn Off.                                                                                                                                                                                                             |
|        | Audible Alarm Button         | Press once to toggle between Audible and Mute Modes.<br>The panel will display the appropriate icon to indicate which mode is enabled.                                                                                                                                              |
|        | Flow Setting Control Buttons | Adjust the oxygen flow setting by pressing the "+" or "-" buttons.                                                                                                                                                                                                                  |
|        | Settings Button              | Press the button to toggle between Pulse and Continuous Modes.<br>Press and hold to view device information, including battery temperature and status, molecular sieve temperature, runtime, and status, device temperature and runtime, as well as firmware and hardware versions. |

## User Interface Buttons and Display Symbols-Continued



Home Screen:

| Icon                                                                                | Description                                                                                                                                             |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3.0                                                                                 | Current Flowrate (range is 1.0-3.0 in Continuous Mode).                                                                                                 |
|    | Battery charge level: <ul style="list-style-type: none"> <li>Battery percentage remaining.</li> <li>Battery time remaining (hours: minutes).</li> </ul> |
|   | Device runtime since powered on: (hours: minutes).                                                                                                      |
|  | Mute Mode                                                                                                                                               |
|  | Audible Mode                                                                                                                                            |

The screen will also display the following icons (continued on next page):



## User Interface Buttons and Display Symbols-Continued

| Icon                                                                               | Description                                                                                                                                                                                                  |
|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | Powered by AC or DC only.                                                                                                                                                                                    |
|   | <ol style="list-style-type: none"> <li>1. Powered by battery only (not plugged in and charging).</li> <li>2. Battery level percentage and time remaining (hours: minutes).</li> <li>3. Device On.</li> </ol> |
|   | <ol style="list-style-type: none"> <li>1. Battery is charging.</li> <li>2. Battery level percentage and estimated time to fully charge the battery.</li> <li>3. Device On.</li> </ol>                        |
|   | <ol style="list-style-type: none"> <li>1. Battery is charging.</li> <li>2. Battery level percentage and estimated time to fully charge the battery.</li> <li>3. Device Off.</li> </ol>                       |
|   | The device has detected an active alarm while in Mute Mode.                                                                                                                                                  |
|   | The device has detected an active alarm while in Audible Mode.                                                                                                                                               |
|   | Runtime of device (hours: minutes) since it was powered on (Example shows 2 hours 35 minutes).                                                                                                               |
|  | Multiple alerts have been detected.<br>Display will scroll to display all alerts.                                                                                                                            |

## User Interface Buttons and Display Symbols-Continued



### Alerts in Audible Mode:

The device has detected an active alarm in Audible Mode.

### Alerts in Mute Mode:

The device has detected an active alarm in Mute Mode.

## **Alerts**

### **Adapter plug/unplug:**

An adapter icon displays when adapter is plugged in and disappears when unplugged. If enabled, an audible alarm will sound.

### **Battery plug/unplug:**

A battery icon displays when battery is connected and disappears when disconnected. If enabled, an audible alarm will sound.

| <b>Alarm Item</b>                | <b>Alarm Condition</b>                                          | <b>System Process</b>                                                                                                      | <b>Display of Screen</b>                                    |
|----------------------------------|-----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|
| Battery Exhausted                | Battery cycle > 500<br>or battery health<br>< 50%               | Alarm only                                                                                                                 | Battery Exhausted<br>Replace Battery                        |
| Low Input Voltage                | Adapter input < 17.0 V                                          | Auto-switch to battery<br>until the adapter input<br>> 18 V                                                                | Low Input Voltage<br>Check Adapter                          |
| Absence of Breath                | No breath detected<br>continuously for more<br>than 60 seconds  | Alarm only                                                                                                                 | Absence of Breath<br>Check Cannula                          |
| Oxygen<br>Concentration<br>< 87% | Concentration < 87%<br>continuously for more<br>than 10 minutes | Alarm only                                                                                                                 | Low O <sub>2</sub> : < 87%<br>Contact Provider              |
| Low Battery                      | 5% ≤ RSOC ≤ 20%<br>Without adapter                              | Alarm only                                                                                                                 | Low Battery<br>Charge Now                                   |
| Oxygen<br>Concentration<br>< 50% | Concentration < 50%<br>continuously for more<br>than 10 minutes | If the condition<br>improves, the alarm will<br>clear automatically<br>If the condition persists,<br>contact your provider | Low O <sub>2</sub> : < 50%<br>Contact Provider              |
| Breath Sensor<br>Fail            | Breath sensor failed                                            | If the condition<br>improves, the alarm will<br>clear automatically<br>If the condition persists,<br>contact your provider | Breath Sensor Fail<br>Check if Cannula<br>Kinked or Blocked |
| Oxygen Sensor<br>Fail            | Oxygen sensor failed                                            | If the condition<br>improves, the alarm will<br>clear automatically<br>If the condition persists,<br>contact your provider | Oxygen Sensor Fail<br>Contact Provider                      |

**Alerts-Continued**

| <b>Alarm Item</b>  | <b>Alarm Condition</b>                                            | <b>System Process</b>           | <b>Display of Screen</b>                                |
|--------------------|-------------------------------------------------------------------|---------------------------------|---------------------------------------------------------|
| Gas Delivery Fail  | No delivery detected after injection                              | Alarm only                      | Gas Delivery Fail<br>Contact Provider                   |
| Gas Obstruction    | Pipe or nasal blocked                                             | Alarm only                      | Gas Obstruction<br>Contact Provider                     |
| Tank Pressure Fail | Tank pressure failed                                              | Alarm only                      | Tank Pressure Fail<br>Contact Provider                  |
| Sieve Bed Fail     | Sieve bed failure or invalid                                      | Auto-shut down after 10 seconds | Sieve Bed Fail<br>Contact Provider                      |
| Compressor Fail    | Compressor failed                                                 | Auto-shut down after 10 seconds | Compressor Fail<br>Contact Provider                     |
| Valve Check Fail   | Valve switch failed                                               | Auto-shut down after 10 seconds | Valve Check Fail<br>Contact Provider                    |
| Cooling Fan Fail   | Cooling fan failed                                                | Auto-shut down after 10 seconds | Cooling Fan Fail<br>Contact Provider                    |
| Battery Depleted   | RSOC < 5% without adapter                                         | Auto-shut down after 10 seconds | Battery Depleted<br>Replace Battery/<br>Connect Adapter |
| System Cold        | System temperature < 32°F/0°C                                     | Alarm only                      | System Cold<br>Move to Warmer Place                     |
| Battery Cold       | Battery temperature < 32°F/0°C                                    | Alarm only                      | Battery Cold<br>Move to Warmer Place                    |
| System Hot         | System temperature > 140°F/60°C                                   | Auto-shut down after 10 seconds | System Hot<br>Move to Cooler Place                      |
| Battery Hot        | Battery temperature > 149°F/65°C                                  | Auto-shut down after 10 seconds | Battery Hot<br>Shut Down,<br>Only Use Adapter           |
| Gas Supply Fail    | Flow or concentration below normal after startup                  | Alarm only                      | Gas Supply Fail<br>Contact Provider                     |
| Sys Startup Fail   | Concentration below 87% continuously > 120 s after system startup | Alarm only                      | System Startup Fail<br>Contact Provider                 |
| Power Supply Fail  | System voltage < 10.5 V                                           | Auto-shut down after 10 seconds | Power Supply Fail<br>Contact Provider                   |

## **Power Supply**

### **Standard Lithium-Ion Battery # TOCBY-2**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is powered by a standard lithium-ion battery.


When fully charged, the TOCBY-2 battery provides up to 6.6 hours of operation. Use the AC or DC Power Supply to recharge the battery; full charging takes no more than 4.5 hours.

### **AC Power Supply # TOCAC**

The AC Power Supply is used to power the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator from an AC power source. It automatically adapts to input voltages ranging from 100 V to 240 V (50-60 Hz), making it compatible with most power sources worldwide.

1. Connect C to D Port.
2. Connect B to the Concentrator.
3. Connect A to the nearest AC power outlet.

### **DC Power Supply # P2DCA-1**

The DC Power Supply is designed for use with the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. The DC power input cable connects directly to a vehicle's cigarette lighter or auxiliary DC power outlet. It accepts an input of 11-16 VDC with a fuse rating of 15 A/125 V, and provides an output of 19 V at 6.3 A.

## **Power Supply-Continued**

### **WARNING**

Use of non-approved the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator accessories may cause damage to the device or user and will void the manufacturer's warranty.

## **Accessories**

### **Nasal Cannula # TOCNC**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator must requires a single lumen Nasal Cannula to provide oxygen to the patient.

The following Nasal Cannula recommended for use with the Oxygen Concentrator.

Nasal Cannula Model: # TOCNC.

Nasal Cannula Manufacturer: HSINER CO., LTD.

### **WARNING**

Nasal Cannulas should not be used by more than one person. DO NOT use a cannula that has a length exceeding 25 ft (7.6 m).

### **CAUTION**

When using a long cannula, the flow setting may need to be increased. Increasing the cannula length may reduce the perceived noise during oxygen bolus delivery.

### **CAUTION**

The Nasal Cannula is designed for disposable use.

### **CAUTION**

Select FDA approved Nasal Cannula (Rhythm Healthcare Item #: TOCNC).

### **CAUTION**

It is recommended to replace the Nasal Cannula after more than 24 hours of continuous use.

## **Accessories-Continued**

### **Carry Bag # TOCCB**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Carry Bag can be securely attached to the cart for convenient transport during daily activities. The front pocket is ideal for carrying the AC Power Supply or personal belongings of the patient.

### **Cart # TOCCRT**

Place the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator inside the Carry Bag, then position the Carry Bag on the Cart base. Use the adjustable Velcro straps to securely fasten the Carry Bag to the Cart handle.

## **Accessories-Continued**

### **Humidifier (not included)**

If your provider has prescribed a humidifier, follow the manufacturer's instructions for use. Attach the humidifier to the oxygen outlet port of the oxygen system. Use of humidifiers not specified for this oxygen system may impair performance and void the warranty.

Select the appropriate humidifier according to the size of the oxygen outlet. The inner diameter of the humidifier connecting pipe should be 4-5 mm.

The following humidifier is recommended for use with the Rhythm Healthcare TOC3 Oxygen Concentrator:

Model: # 30506

Registered Establishment Number: 3003862188

Humidifier Manufacturer: HSINER CO., LTD.

### **Setting up the Humidifier:**

1. Add distilled water to the humidifier. See Figure 1.
2. Secure the humidifier to the device with the strap. See Figure 2.
3. Connect the TOC3 Oxygen Outlet tubing to the Humidifier Oxygen Inlet with the humidifier. See Figure 3.
4. Connect the Nasal Cannula to the Humidifier outlet. See Figure 4.

**Figure 1**

**Figure 2**

**Figure 3**

TOC3 Oxygen  
Outlet

Humidifier:  
Oxygen Inlet

**Figure 4**

### **WARNING**

DO NOT use the Humidifier in Pulse Mode.

### **CAUTION**

DO NOT overfill the Humidifier, if used. Fill only with water up to the level indicated by the manufacturer.

### **CAUTION**

The Nasal Cannula must be attached when used with the Humidifier.

## **The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Accessories**

| <b>Item Number</b> | <b>Description</b>                |
|--------------------|-----------------------------------|
| TOC3               | Transportable Oxygen Concentrator |
| TOCNC              | Nasal Cannula                     |
| TOCCB              | Carry Bag                         |
| TOCAC              | AC Power Supply                   |
| P2DCA-1            | DC Power Supply                   |
| TOCIF-1            | Small Intake Filter               |
| TOCIF-2            | Large Intake Filter               |
| TOCBY-2            | Battery                           |
| TOCCRT             | Cart                              |

## **General Operation**

- 1. Find a well ventilated location to place the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Ensure that the unit is turned off.**

Ensure the Intake and exhaust has clear access. The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator should be placed where alarms can be heard.

### **WARNING**

Do not use the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator in the presence of flammable anesthetics, detergents, or other chemical vapors.

### **CAUTION**

The oxygen concentrator requires 30 minutes to warm up from the minimum storage temperature before it is ready for use. Similarly, it requires 30 minutes to cool down from the maximum storage temperature before it is ready for use.

### **CAUTION**

Do not block the Air Intake (Pre-Filter Cover) or Air Exhaust when operating the equipment. Blocked air circulation or proximity to the heat source can cause internal heat build-up, shut down, or damage to the concentrator.

### **CAUTION**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator can be used 24 hours a day while traveling. Regular use is recommended to preserve optimal performance and extend the life of the sieve beds.

### **CAUTION**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is shipped from the factory with the battery removed and not installed.

### **CAUTION**


Do not operate the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator without the Small and Large Intake Filters in place. Operating the device without these filters in place can damage the device.

## **General Operation-Continued**

### **2. Install the battery.**

Slide the battery into place until the latch clicks into place. An audible click will confirm that the battery is securely installed.

### **3. Connect the AC or DC Power Supply to the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.**

#### **CAUTION**

Only use the supplied power cord with the power port. Do not insert any other items into the power port. Do not use extension cords.

#### **CAUTION**

Power supply is not waterproof. Do not disassemble the power supply.


#### **CAUTION**

When unplugging the power from the wall outlet, also disconnect it from the concentrator to avoid unnecessary battery drain.

## **General Operation-Continued**

### **4. Connect the Nasal Cannula to the Nozzle Fitting.**

Nozzle fitting is located on the top of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Connect the Nasal Cannula to the nozzle fitting on the device (pictured).

#### **CAUTION**

Ensure that the cannula is routed to prevent it from being pinched or kinked to avoid a disruption of oxygen flow.

#### **CAUTION**

The proper placement and positioning of the prongs of the Nasal Cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.

### **5. Press the On/Off button to turn on the concentrator.**

At startup, the device will beep and the indicator light will flash.

"Welcome" will appear on the display at start up. A five-minute warm-up period will begin, during which the oxygen concentration will gradually increase but may not immediately reach the specified level.

In certain conditions, such as exposure to very cold temperatures during storage or use, a longer warm-up time may be necessary.

#### **CAUTION**

Oxygen concentration may not reach specification during the five-minute warm up time.

#### **CAUTION**


30 seconds after startup, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will enter Auto-Pulse Mode. Inhalation will not be detected or supported during this initial 30 second period.

#### **CAUTION**

Ensure the AC Power Supply is in a well-ventilated area. During operation, the power supply may get hot. Ensure the power supply is cool before handling.

## General Operation-Continued

6. Wear the Nasal Cannula on your face and breathe through your nose.


7. Select Pulse or Continuous Mode and adjust flow to the prescribed level.

Press the Settings button to switch between Pulse and Continuous Modes.

Press the "+" or "-" buttons to adjust to your prescribed flow rate.

## **General Operation-Continued**

### **In Pulse Mode:**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is designed to detect when you are breathing through the cannula. If you haven't started breathing from it yet, the device will automatically deliver a pulse of oxygen approximately every three seconds. Once you begin breathing through the cannula, the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will synchronize with your breath, delivering oxygen pulses in response to your inhalation. As your breathing rate changes, the device senses these changes and adjusts the pulse volume accordingly for your next breath.

### **WARNING**

**DO NOT use a humidifier in Pulse Mode.**

### **In Continuous Mode:**

A constant flow rate of oxygen will be provided.

### **WARNING**

If you feel discomfort using the device, consult your provider immediately.

### **CAUTION**

If the oxygen level drops below recommended levels, a Low O<sub>2</sub>: <87% alert will display on the screen. If the alarm persists, contact your provider.

### **CAUTION**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator will notify you with an audible alert (if enabled) and display "Absence Of Breath" if no breath has been detected for 60 seconds. After 60 seconds, the device will enter Auto Pulse Mode until breath is detected. Once breath is detected, the device will resume normal delivery of oxygen.

### **Note:**

To cut off power, unplug the input cord from the AC Power Supply (such as a wall outlet or DC Power Supply), then disconnect it from the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

## **Troubleshooting**

The table below lists some common problems and solutions. If you can't resolve a problem, please contact your provider.

| <b>Problem</b>                   | <b>Possible Cause</b>                      | <b>Recommended Solution</b>                                                                                                                                           |
|----------------------------------|--------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Device Won't Turn On             | Battery is not installed correctly.        | Remove the battery and re-install it correctly.                                                                                                                       |
|                                  | Battery is depleted.                       | Use the AC or DC Power Supply to operate the device (with the battery inserted) to recharge the battery. If this does not resolve the problem, contact your provider. |
|                                  | AC Power Supply is not connected properly. | Check AC power connection and verify solid green light on power supply.                                                                                               |
|                                  | DC Power Supply is not connected properly. | Check DC connection and at cigarette lighter or auxiliary DC power source.                                                                                            |
| No Oxygen                        | The device is not turned on.               | Turn on the concentrator.                                                                                                                                             |
|                                  | Cannula is kinked or obstructed.           | Check cannula and its connection to the oxygen outlet port.                                                                                                           |
|                                  | Equipment failure.                         | Contact your provider.                                                                                                                                                |
| Oxygen Not at Full Concentration | The device is warming up.                  | Wait 5 minutes for the device to warm-up. If the problem is not solved, please contact your provider.                                                                 |
|                                  | The sieve beds may require servicing.      | Contact your provider to change the sieve beds.                                                                                                                       |
| Alarm Occurs                     | Refer to previous section-Alerts.          | Refer to previous section-Alerts.                                                                                                                                     |

## **Cleaning the Case**

The outside case should be cleaned using a damp cloth with a solution of mild detergent and water when necessary.

### **CAUTION**

Do not allow liquids into any of the controls, the interior of the case, or the oxygen tubing connector. If this occurs, contact your provider for assistance.

### **WARNING**

Do not use alcohol, isopropyl alcohol, ethylene chloride or petroleum-based cleaners on the cases or on the filters.

## **Cannula Replacement**

The Nasal Cannula is disposable. Contact your provider should you need Nasal Cannula replacements. Follow the cannula manufacturer's instructions.

### **CAUTION**

Nasal Cannula should be FDA approved and the length should be less than 25 ft (7.6 m).

### **CAUTION**


Nasal Cannula should be rated for 5 liters per minute to ensure proper patient usage and oxygen delivery.

## **Filter Cleaning and Replacement**

### **Pre-Filter Cover # TOCSP-06**

The Pre-Filter Cover is designed for adequate air flow through the device and to ensure clean air is delivered into the compressor.

The Pre-Filter Cover must be cleaned weekly to ensure adequate air flow. Clean with a mild liquid detergent and water. Ensure that the Pre-Filter Cover is completely dry before use.


### **CAUTION**

It may be necessary to clean the Pre-Filter Cover more often in dusty or polluted environments/conditions.

### **Intake Filters**


The filters are designed to allow proper airflow through the bottom of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator includes two intake filters. Clean the filters with a mild liquid detergent and water. Ensure they are completely dry before reinstalling and using them.


### **Replacing the Small Intake Filter (TOCIF-1)**

1. Remove the Filter Box from the bottom of the device.

## **Filter Cleaning and Replacement-Continued**

2. Remove the Small Intake Filter from the Intake Chamber.

3. Insert a new Small Intake Filter into the chamber.
4. Reinstall the Filter Box.

### **CAUTION**

Small Intake Filter (TOCIF-1) should be replaced every 6 months and cleaned once a week. If the device is used in a dusty environment, the Small Intake Filter may need to be replaced more frequently. Do not operate the device without both intake filters installed.

### **Replacing the Large Intake Filter (TOCIF-2)**

1. Remove the Filter Box from the bottom of the device.


2. Remove the Filter Plate.

## **Filter Cleaning and Replacement-Continued**

3. Remove the Large Intake Filter from the Filter Box.

4. Insert a new Large Intake Filter.
5. Place the Filter Plate back in the Filter Box.
6. Reinstall the Filter Box.

### **CAUTION**

The Large Intake Filter (TOCIF-2) should be replaced every year and cleaned once a month. If the device is used in a dusty environment, the filter may need to be replaced more frequently. Do not operate the device without intake filters installed.

**Contact your provider for replacement Pre-Filter Covers and Intake Filters.**

## **Battery Care and Maintenance**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator Lithium-Ion Battery requires special care to ensure proper performance and long life. Only use the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator # TOCBY-2 battery with your concentrator.

### **CAUTION**

Keep liquids away from battery. If battery gets wet, please stop using immediately and dispose of battery properly.

### Battery Replacement

1. Press down on the latch and slide the battery out.
2. Slide battery into place until the latch clicks into place.

### Effect of Temperature on Battery Performance

To maximize battery run-time, operate the device in temperatures between 41°F and 95°F (5°C and 35°C). Battery lifespan and performance are strongly influenced by the temperature at which the battery is charged.

### CAUTION

The room temperature should not exceed 75°F (24°C) when batteries are being charged.

### Battery Time Remaining

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator continuously shows the estimated remaining battery time. Please note that this is only an approximation, and the actual remaining time may differ.

### CAUTION

Store the battery in a cool, dry location. To preserve its health and extend its lifespan, charge the battery to 80% every three months if the device will be unused for an extended period. Also, be sure to remove the battery from the device during storage.

### BATTERIES SHOULD NOT BE LEFT DORMANT FOR MORE THAN 90 DAYS AT A TIME.

### CAUTION

If the device is not used for an extended period of time, please remove the battery from the device and charge it every three months.

## **Disposal of Equipment and Accessories**

Dispose of and recycle Rhythm Healthcare TOC3 Transportable Oxygen Concentrator accessories in accordance with local regulations. The battery contains lithium-ion cells and must be recycled properly-do not incinerate.

## **Replacement Items**

| <b>Part number</b> | <b>Description</b>  |
|--------------------|---------------------|
| TOCBY-2            | Battery             |
| TOCIF-1            | Small Intake Filter |
| TOCIF-2            | Large Intake Filter |

If you need assistance, please contact your provider.

**Specifications**

|                                               |                                                                                                                                                                     |                    |    |    |    |     |     |     |     |     |     |  |
|-----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|----|----|----|-----|-----|-----|-----|-----|-----|--|
| Mode                                          | Setting                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
| Continuous Flow                               | 1.0 LPM                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
|                                               | 1.5 LPM                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
|                                               | 2.0 LPM                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
|                                               | 2.5 LPM                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
|                                               | 3.0 LPM                                                                                                                                                             |                    |    |    |    |     |     |     |     |     |     |  |
|                                               | Deliverable flow accuracy:<br>The flow error is less than 10% or 0.2 L/min, whichever is greater.                                                                   |                    |    |    |    |     |     |     |     |     |     |  |
| Mode                                          | Setting                                                                                                                                                             | Flowrate L/Min     |    |    |    |     |     |     |     |     |     |  |
| Pulse Dose                                    | 1                                                                                                                                                                   | 0.2                |    |    |    |     |     |     |     |     |     |  |
|                                               | 2                                                                                                                                                                   | 0.4                |    |    |    |     |     |     |     |     |     |  |
|                                               | 3                                                                                                                                                                   | 0.6                |    |    |    |     |     |     |     |     |     |  |
|                                               | 4                                                                                                                                                                   | 0.8                |    |    |    |     |     |     |     |     |     |  |
|                                               | 5                                                                                                                                                                   | 1.0                |    |    |    |     |     |     |     |     |     |  |
|                                               | 6                                                                                                                                                                   | 1.2                |    |    |    |     |     |     |     |     |     |  |
|                                               | 7                                                                                                                                                                   | 1.4                |    |    |    |     |     |     |     |     |     |  |
|                                               | 8                                                                                                                                                                   | 1.6                |    |    |    |     |     |     |     |     |     |  |
|                                               | 9                                                                                                                                                                   | 1.8                |    |    |    |     |     |     |     |     |     |  |
|                                               | 10                                                                                                                                                                  | 2.0                |    |    |    |     |     |     |     |     |     |  |
| Flow Control<br>Settings and Pulse<br>Volumes | Setting                                                                                                                                                             | 1                  | 2  | 3  | 4  | 5   | 6   | 7   | 8   | 9   | 10  |  |
|                                               | Breath<br>Rate                                                                                                                                                      | Pulse Volumes (mL) |    |    |    |     |     |     |     |     |     |  |
|                                               | 10                                                                                                                                                                  | 20                 | 40 | 60 | 80 | 100 | 120 | 140 | 160 | 180 | 200 |  |
|                                               | 15                                                                                                                                                                  | 13                 | 27 | 40 | 53 | 67  | 80  | 93  | 107 | 120 | 133 |  |
|                                               | 20                                                                                                                                                                  | 10                 | 20 | 30 | 40 | 50  | 60  | 70  | 80  | 90  | 100 |  |
|                                               | 25                                                                                                                                                                  | 8                  | 16 | 24 | 32 | 40  | 48  | 56  | 64  | 72  | 80  |  |
|                                               | 30                                                                                                                                                                  | 7                  | 13 | 20 | 27 | 33  | 40  | 47  | 53  | 60  | 67  |  |
|                                               | 35                                                                                                                                                                  | 6                  | 11 | 17 | 23 | 29  | 34  | 40  | 46  | 51  | 57  |  |
|                                               | 40                                                                                                                                                                  | 5                  | 10 | 15 | 20 | 25  | 30  | 35  | 40  | 45  | 50  |  |
|                                               | Deliverable flow accuracy: $\pm 15\%$ at STPD*<br>$\pm 25\%$ over the rated environmental range<br>*STPD is 101.3 kPa at an operating temperature of 68°F/20°C, dry |                    |    |    |    |     |     |     |     |     |     |  |
| Dimensions<br>(Length x Width x Height)       | 11.1" L x 6.7" W x 15" H<br>(28.2 cm L x 17.0 cm W x 38.1 cm H)                                                                                                     |                    |    |    |    |     |     |     |     |     |     |  |
| Weight                                        | 17.2 lbs (7.8 kg) (With TOCBY-2)                                                                                                                                    |                    |    |    |    |     |     |     |     |     |     |  |

**Specifications-Continued**

|                                                                   |                                                                                                                                                                                         |                                                             |
|-------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|
| User Interface                                                    | 2.8" LCD color display screen                                                                                                                                                           |                                                             |
| Sound Level                                                       | 42 dBA at Pulse Setting 2<br>45 dBA at Continuous Setting 2 LPM (MDS-Hi)* <sup>1, 2</sup>                                                                                               |                                                             |
| Audio Alarm                                                       | ≥45 dBA                                                                                                                                                                                 |                                                             |
| Time from Turning on the Concentrator to Reach Stated Performance | 5 minutes                                                                                                                                                                               |                                                             |
| Oxygen Concentration                                              | 87%-96% (all settings)                                                                                                                                                                  |                                                             |
| Breathing Frequency                                               | 10 to 40 bpm                                                                                                                                                                            |                                                             |
| Inspiratory Trigger Sensitivity                                   | ≤-4.72" H <sub>2</sub> O                                                                                                                                                                |                                                             |
| Time Delay from Onset of Inspiration                              | <10 ms                                                                                                                                                                                  |                                                             |
| Delivery Pressure at The Device Outlet                            | Pulse Dose: 20 PSI (+10%)<br>Continuous Flow: 5 PSI (+10%)                                                                                                                              |                                                             |
| Power                                                             | AC Power Supply                                                                                                                                                                         | AC Input: 100 to 240 VAC 50 to 60 Hz                        |
|                                                                   | DC Power Supply                                                                                                                                                                         | DC Input: 11-16 VDC<br>DC Output: 19 V = 6.3 A              |
|                                                                   | Rechargeable Battery                                                                                                                                                                    | TOCBY-2: Voltage: 14.4 VDC<br>Rated capacity: 6.7 Ah+6.7 Ah |
| Battery Charging Time                                             | TOCBY-2: Up to 4.5 hours<br>3.4 Hours at Device Standby Mode when over 90% fully charged by AC power supply.<br>4.5 Hours at Device Standby Mode when fully charged by AC power supply. |                                                             |
| Environmental Ranges Intended for Operation                       | Temperature: 41°F to 104°F (5°C to 40°C)<br>Humidity: 10% to 90% Non-condensing<br>Altitude: 0 to 10,000 ft (0 to 3,048 meters, 70 kPa to 106 kPa)                                      |                                                             |
| Environmental Ranges Intended for Shipping and Storage            | Temperature: -4°F to 158°F (-20°C to 70°C)<br>Humidity: 5% to 90% Non-condensing<br>Store in a dry environment<br>Altitude: 0 to 10,000 ft (0 to 3,048 meters, 70 kPa to 106 kPa)       |                                                             |

\* 1 According to ISO 80601-2-69 the Maximum system sound pressure level of the device is 49 dBA at Pulse Setting 2 and 55 dBA at Continuous Setting 3 LPM.

\* 2 Data according to test method 14-1 03/2007 MDS-Hi.

**Specifications-Continued**

|                  |                                                                                                    |                 |                                                                                                                                                                                                       |
|------------------|----------------------------------------------------------------------------------------------------|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Battery Duration | TOCBY-2                                                                                            | Continuous Flow | 1.0 LPM: 3.1 h<br>1.5 LPM: 2.5 h<br>2.0 LPM: 1.9 h<br>2.5 LPM: 1.4 h<br>3.0 LPM: 1.2 h                                                                                                                |
|                  |                                                                                                    | Pulse Dose      | Setting 1: 6.6 h<br>Setting 2: 4.9 h<br>Setting 3: 3.9 h<br>Setting 4: 3.6 h<br>Setting 5: 3.3 h<br>Setting 6: 3.0 h<br>Setting 7: 2.6 h<br>Setting 8: 2.3 h<br>Setting 9: 2.0 h<br>Setting 10: 1.8 h |
| Power Modes      | Battery: No Limitation                                                                             |                 |                                                                                                                                                                                                       |
|                  | AC Power Supply: No Limitation                                                                     |                 |                                                                                                                                                                                                       |
|                  | DC Power Supply: Limited to Settings 1-8 in Pulse Mode<br>Limited to 1- 1.5 LPM in Continuous Mode |                 |                                                                                                                                                                                                       |
| Operating Time   | 24-hours a day when connected to an external AC or DC power source.                                |                 |                                                                                                                                                                                                       |
| Transportation   | Keep Dry, Handle With Care.                                                                        |                 |                                                                                                                                                                                                       |

**Classifications**

|                                                                          |                                                                                                               |
|--------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Mode of Operation                                                        | Continuous Duty                                                                                               |
| Type of Protection Against Electrical Shock                              | Class II                                                                                                      |
| Degree of Protection to Concentrator Components Against Electrical Shock | Type BF<br>Not intended for cardiac application                                                               |
| Degree of Protection to Concentrator Components Against Ingress of Water | IP22-Protects against solid objects over 12.5 mm and direct sprays of water up to 15° of vertical (IEC 60529) |

## **Standards Compliance**

The device is designed to conform to the following standards:

- IEC 60601-1: 2005/ADM1: 2012/ADM2: 2020, Medical electrical equipment-Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2: 2014/AMD1: 2020, Medical electrical equipment-Part 1-2: General requirements for basic safety and essential performance-Collateral Standard: Electromagnetic disturbances-Requirements and tests
- IEC 60601-1-6: 2010/AMD1: 2013/AMD2: 2020, Medical electrical equipment-Part 1-6: General requirements for basic safety and essential performance-Collateral standard: Usability
- IEC 60601-1-8: 2006/AMD1: 2012/AMD2: 2020, Medical electrical equipment-Part 1-8: General requirements for basic safety and essential performance-Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
- IEC 60601-1-11: 2015/AMD1: 2020, Medical electrical equipment-Part 1-11: General requirements for basic safety and essential performance-Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- ISO 80601-2-67: 2020, Medical electrical equipment-Part 2-67: Particular requirements for basic safety and essential performance of oxygen-conserving equipment
- ISO 80601-2-69: 2020, Medical electrical equipment-Part 2-69: Particular requirements for the basic safety and essential performance of oxygen-concentrator equipment
- ISO 18562-1: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications-Part 1: Evaluation and testing within a risk management process
- ISO 18562-2: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications-Part 2: Tests for emissions of particulate matter
- ISO 18562-3: 2017, Biocompatibility evaluation of breathing gas pathways in healthcare applications-Part 3: Tests for emissions of volatile organic compounds (VOCs)
- ISO 10993-1: 2018, Biological evaluation of medical devices-Part 1: Evaluation and testing within a risk management process
- ISO 10993-5: 2009, Biological evaluation of medical devices-Part 5: Tests for in vitro cytotoxicity
- ISO 10993-10: 2021, Biological evaluation of medical devices-Part 10: Tests for skin sensitization
- ISO 10993-23: 2021, Biological evaluation of medical devices-Part 23: Tests for irritation

## **EMC Information**

The device has been designed to meet EMC standards throughout its Service Life.

Guidance and Manufacturer's Declaration-Electromagnetic Immunity:

The Concentrator is intended for use in the electromagnetic environment specified below.

The users of the Concentrator should make sure that it is used in such an environment.

| <b>Immunity Test</b>                                                         | <b>IEC 60601 Test Level</b>                                                                                                                                           | <b>Compliance Level</b>                                                                                                                                               | <b>Electromagnetic Environment-Guidance</b>                                                                                                                                                                                                                          |
|------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic Discharge (ESD)<br>IEC 61000-4-2                               | ±8 kV Contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV Air                                                                                                                   | ±8 kV Contact<br>±2 kV, ±4 kV, ±8 kV,<br>±15 kV Air                                                                                                                   | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                        |
| Electrical Fast Transient/Burst IEC 61000-4-4                                | ±2 kV for Power Supply Lines<br>±1 kV for Input/Output Lines                                                                                                          | ±2 kV for Power Supply Lines<br>±1 kV for Input/Output Lines                                                                                                          | Mains power quality should be that of a typical home or hospital environment.                                                                                                                                                                                        |
| Surge<br>IEC 61000-4-5                                                       | ±1 kV Line to Line                                                                                                                                                    | ±1 kV Line to Line                                                                                                                                                    | Main power quality should be that of a typical home or hospital environment.                                                                                                                                                                                         |
| Voltage Dips, Short Interruptions on Power Supply Input Lines IEC 61000-4-11 | <5% $U_T$ (>95% Dip in $U_T$ ) for 0.5 Cycle at 45 degree increments<br>70% $U_T$ (30% Dip in $U_T$ ) for 0.5 Seconds<br><5% $U_T$ (>95% Dip in $U_T$ ) for 5 Seconds | <5% $U_T$ (>95% Dip in $U_T$ ) for 0.5 Cycle at 45 degree increments<br>70% $U_T$ (30% Dip in $U_T$ ) for 0.5 Seconds<br><5% $U_T$ (>95% Dip in $U_T$ ) for 5 Seconds | Mains power quality should be that of a typical home or hospital environment. If the users of the device be required continued operation during power mains interruptions, it is recommended that device be powered from an uninterruptible power supply or battery. |
| Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8                      | 30 A/m                                                                                                                                                                | 30 A/m                                                                                                                                                                | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical home or hospital environment.                                                                                                                                  |

Note:  $U_T$  is the A.C. mains voltage prior to application of the test level.

**EMC Information-Continued**

|                               |                                                                |                                                                |                                                                                                                                                                                                                                                                                                                                                          |
|-------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz                                    | 3 Vrms<br>150 kHz to 80 MHz                                    | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended 12 inches separation distance.<br><br>Interference may occur in the vicinity of equipment marked with the following symbol:  |
| Radiated RF<br>IEC 61000-4-3  | 6 Vrms<br>Amateur Radio & ISM Bands between 150 kHz and 80 MHz | 6 Vrms<br>Amateur Radio & ISM Bands between 150 kHz and 80 MHz |                                                                                                                                                                                                                                                                                                                                                          |
| Radiated RF<br>IEC 61000-4-3  | 10 V/m<br>80 MHz to 2.7G Hz                                    | 10 V/m                                                         |                                                                                                                                                                                                                                                                                                                                                          |

**CAUTION**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator is may be used in home, institution and various mobile environments.

**CAUTION**

Electromagnetic interference could disturb performance of the device (oxygen concentration and flow), resulting in an alarm.

When used in the above environment, it is also necessary to avoid exposure to known sources of Electro-Magnetic Interference, including MRI, RFID, diathermy, and electrocautery.

**CAUTION**

Use of accessories (Battery, AC Power Supply, DC Power Supply) other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic, immunity of this equipment and result in improper operation.

**CAUTION**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches to any part of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator, including cables supplied by the manufacturer. The use of this device closer than the recommended distance can result in a degradation of the performance of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator.

## **EMC Information-Continued**

### **CAUTION**

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

## **Guidance and Manufacturer's Declaration-Electromagnetic Emissions**

The Concentrator is intended for use in the electromagnetic environment specified below. The user of the Concentrator should make sure it is used in such an environment.

| <b>Emissions test</b>                                | <b>Compliance</b> | <b>Electromagnetic Environment-Guidance</b>                                                                                                                                                                              |
|------------------------------------------------------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF Emissions CISPR 11                                | Group 1           | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                          |
| RF Emissions CISPR 11                                | Class B           | The device is suitable for uses in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic Emissions IEC 61000-3-2                     | Class A           |                                                                                                                                                                                                                          |
| Voltage Fluctuations/Flicker Emissions IEC 61000-3-3 | Complies          |                                                                                                                                                                                                                          |

**Essential Performance List**

| <b>Essential Performances</b>                                  | <b>Relevant Requirements or Conditions</b>                                                             |
|----------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
| Accuracy of concentration                                      | Within 87% - 96% (at all settings)                                                                     |
| Accuracy of oxygen delivery                                    | Within 25% as indicated in the instruction for use                                                     |
| Power Supply Failure Alert                                     | System voltage < 10.5 V                                                                                |
| Internal electrical power source nears depletion               | RSOC < 5% without adapter                                                                              |
| Low Oxygen Concentration Alert                                 | The oxygen concentration remains below 87% continuously for more than 600 seconds                      |
| Start Up Conditions Alert                                      | The oxygen concentration remains below 87% continuously for more than 120 seconds after system startup |
| Triggering signal absence alarm condition (Only in Pulse Mode) | No breath has been detected continuously for more than 60 seconds.                                     |
| Gas Supply Failure Alert                                       | Flow or concentration below normal after injection                                                     |

## **Warranty**

The Rhythm Healthcare TOC3 Transportable Oxygen Concentrator warranty covers the repair or replacement of the unit. The warranty term is 36 months from the date of shipment. Please contact us by telephone or email to return defective equipment under warranty and to resolve any problems. Our trained technicians will help you with any questions or problems with your the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. Please make sure that your returned equipment is packaged safely for transportation, if possible, in its original packaging to avoid damages during shipping. Excluded from the warranty are damages caused by improper usage. Also excluded are replacements of batteries, disposable parts, and consumables. Sieve bed, filters, batteries are expressly excluded from the 36 months warranty, except as provided below:

| Description                                                           | Period      |
|-----------------------------------------------------------------------|-------------|
| Oxygen Concentrator                                                   | 3 Years     |
| Accessories (Battery, Carry Bag, AC Power Supply and DC Power Supply) | 1 Year      |
| Disposables (Nasal Cannula, Filters)                                  | No Warranty |
| Sieve Bed                                                             | 2 Years     |

Further damage compensation claims of any kind, particularly owing to breach of obligations and unpermitted handling, as well as claims on repayment of expenses paid in vain, are not included in the warranty; the same shall apply to claims on repayments of consequential harm caused by a defect.

Any further claims are excluded in this warranty. The aforementioned limitations do not apply to claims on damages from harm to life, body, or health or attributable to intent or gross negligence, or the product liability law.

This warranty does not cover damage or injury whether to the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator or to personal property or persons caused by accident, misuse, abuse, negligence, failure to install in accordance with Rhythm's installation instructions, failure to operate under conditions of normal use and in accordance with the terms of the operating manual and instructions, failure to maintain in accordance with the applicable service manuals, or alteration or any defects not related to materials or workmanship of the Rhythm Healthcare TOC3 Transportable Oxygen Concentrator. This warranty does not cover damage which may occur in shipment. This warranty does not apply to any product or individual part of a product that may have been repaired or altered by anyone other than Rhythm Healthcare or an authorized Rhythm Healthcare service center. This warranty does not apply to any product which is not purchased new.

### **Adverse Reactions**

MedWatch is the Food and Drug Administration's (FDA) program for reporting serious reactions, product quality problems, therapeutic inequivalence/failure, and product use errors with human medical products, including drugs, biologic products, medical devices, dietary supplements, infant formula, and cosmetics.

If you think you or someone in your family has experienced a serious reaction to a medical product, you are encouraged to take the reporting form to your doctor. Your health care provider can provide clinical information based on your medical record that can help FDA evaluate your report.

However, we understand that for a variety of reasons, you may not wish to have the form filled out by your health care provider, or your health care provider may choose not to complete the form. Your health care provider is not required to report to the FDA. In these situations, you may complete the Online Reporting Form yourself.

You will receive an acknowledgement from FDA when your report is received. Reports are reviewed by FDA staff. You will be personally contacted only if we need additional information.

### **Submitting Adverse Event Reports to FDA**

Use one of the methods below to submit voluntary adverse event reports to the FDA at <http://www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home>

Consumer Reporting Form FDA 3500B.

Follow the instructions on the form to either fax or mail it in for submission. For help filling out the form, see MedWatchLearn.

The form is available at <http://www.fda.gov/downloads/aboutFDA/reportsmanualsforms/forms/ucm349464.pdf>.

Call FDA at 1-800-FDA-1088 to report by telephone. Reporting Form FDA 3500 is commonly used by health professionals.

The form is available at <http://www.fda.gov/downloads/aboutFDA/reportmanualsforms/forms/ucm163919.pdf>.

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Rhythm Healthcare, LLC