




LR200SF
LR300SF
LR400SF
LR500SF
LR600SF

SINGLE FOLD ALUMINUM RAMPS

600 LBS
WEIGHT CAPACITY



USER INSTRUCTIONS



Determine the recommended degree of incline for your assistive device (scooter, wheelchair, etc) by referring to your owner's manual

- Measure the distance from the top step or landing to the ground (RISE).
 - Refer to Incline Chart to find the proper ramp length.
 - Always use a qualified assistant when using a ramp.
1. Your ramp has a 600 pound weight capacity. Please remember this is determined by adding your weight to that of your assistive device.
 2. Always inspect the surface of the ramp for anything that could be slippery such as water, ice, oil, etc. Clean surface if anything is found. Also, remove any debris such as branches, twigs or rubbish before use. While inspecting ramp, always check for worn or damaged parts. If any problems are found, do NOT use the ramp. Contact your dealer for assistance.
 3. Always secure the top flange of your ramp on a level surface and make sure the top edge is overlapping the step by at least 1 1/2 inches. Make sure the bottom of the ramp rests securely on the ground.

4. Always secure the pins firmly in the holes at the top of the ramp. Holes will need to be drilled into the surface in which you are using the ramp to properly secure it and maintain your safety.
5. Do not exceed the recommended manufacturer's degree of incline when using your ramp. Doing so may cause your power chair or scooter to tip forwards or backwards resulting in harm to the patient.
6. Do not add or modify your ramp in any way. THIS WILL VOID YOUR WARRANTY.
7. Frequently inspect your ramp to make sure that the bolts, nuts, and hinges are tight and properly secured.
8. When using your ramp for access into a van: NEVER ride your assistive device into or out of the vehicle! Manually push the assistive device into and out of the vehicle.
9. Always use the single-fold ramp in the full open position only!
10. When using the ramp with a wheelchair always have a qualified attendant with you and ALWAYS go down the ramp backwards, it is highly recommended that a lap belt be worn for safety.
11. Please check the ramp usage chart, to determine which ramp length is best in your application. Failure to do this may cause injury.

WARRANTY

Rhythm Healthcare, LLC ramps are covered with a three (3) year transferable warranty. Please see warranty certificate for details.

Rhythm Healthcare
Clearwater, FL 33762

contactus@rhythmhc.com
Website: www.rhythmhc.com
(877) 843-6464