


Exploder Drawing



NO.	item Number	Part Name	QTY
1	SEHCFT	HeadBoard Complete for the Semi Electric Bed, 1 ea	1
2	FEHCWM	HeadSpring Complete without Motor for the Full Electric Low Be & the Semi Electric Bed, 1 ea	1
3	FEFCWM	FootSpring Complete without Motor for the Full Electric Low Bed & the Semi Electric Bed, 1 ea	1
4	SEFCFT	FootBoard Complete for the Semi Electric Bed, 1 ea	1
5	FEGOFF	Gearbox of Foot for the Full Electric Low Bed & the Semi Electric Bed, 1 ea	1
6	SEHCRA	Hand Crank for the Semi Electric Bed, 1 ea	1
7	FECWLF	Caster with Locking for the Full Electric Low Bed & the Semi Electric Bed, 1 ea	2
8	SEDSFT	Drive Shaft for the Semi Electric Bed, 1 ea	1
9	SEMFTS	Motor for the Semi Electric Bed, 1 ea	2
10	SEHCON	Hand Controller for the Semi Electric Bed, 1 ea	1

11	FECWOLF	Caster without Locking for the Full Electric Low Bed &the Semi Electric Bed, 1 ea	2
12	FEGOHF	Gearbox of Head for the Full Electric Low Bed & the Semi Electric Bed, 1 ea	1

12



11