

# POC BATTERY CHARGER



ITEM #	DESCRIPTION
P2PBC090	Battery Charger

## FEATURES

- Fast charging as 2.5-4 hours.
- Easy operation.
- Conveniently charge the battery while the POC device is in use (Compatible with for the P2, P2-E6, P2-E7, P2-E)

Green indicator lights flashing when charging, Five green solid lights indicate the battery is fully charged.

## WARRANTY

Limited Lifetime Warranty is provided against any factory defects in material or workmanship. One (1) year warranty against defective materials or workmanship on all other components, assemblies and parts. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and/or unauthorized modifications.

For questions, please email us at [Contactus@rhythmhc.com](mailto:Contactus@rhythmhc.com). Visit [www.rhythmhc.com](http://www.rhythmhc.com) for more information about our products and accessories.