

COMBO WHEELCHAIR/ TRANSPORT CHAIR

300 lb Weight Capacity



ITEM #	DESCRIPTION
L3418SLN	18" Silver Combo Wheelchair/Transport Chair
L3420SLN	20" Silver Combo Wheelchair/Transport Chair

FEATURES

- 18"W x 16"D Seat (18" Chair) and 20"W x 16"D (20" Chair)
- Standard Self-Propelled Wheelchair And a Transport Chair All In One
- Quick-release 24" wheels
- Aluminum Frame
- Aluminum crossbar
- Detachable Armrest With New Armrest Pad
- Two Sets of Wheel Locks
- Back Folds Down With Deluxe Back Release
- Nylon Upholstery And Arm Pad
- Swing-away Footrests With Composite Foot Plates and Heel Loops
- 8" Casters
- Aluminum Brake
- Includes Anti-tippers
- Tool-free

SPECIFICATIONS

Back height: 18.1"

Overall Height: 39"

Overall Length w/ Riggings: 41.3"

Seat To Floor Height: 18.9"

Width Open/Closed: 29.1"/14.2"

Product weight: 37.5 lbs

V1: 1/12/23

COMBO WHEELCHAIR/ TRANSPORT CHAIR

300 lb Weight Capacity



LIMITED LIFETIME WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.