



LM5BA 5 LPM STATIONARY OXYGEN CONCENTRATOR

ITEM #	DESCRIPTION
LM5BA	5 LPM Stationary Oxygen Concentrator

FEATURES

- User friendly control panel
- Alarms for temperature, oxygen purity, high and low pressure, power failure, and blockage
- Quick & easy maintenance
- Low noise level
- High output pressure

SPECIFICATIONS

Power Supply	AC115v/60Hz
Power Consumption	<365W
Output Pressure	40-60kPa
Oxygen Flow	5LPM
Oxygen Purity	87%-96%
Equip Classification	Class II Type B
Unit Dimensions	16.9"W x 14.4"L x 25.6"H
Product Weight	35.25 lbs

WARRANTY

Three Year Limited Warranty

V2: 10/12/23