

## P2-E7 PORTABLE OXYGEN CONCENTRATOR



| ITEM # | DESCRIPTION                            |  |
|--------|--|--|
| P2-E7  | 7 Setting Portable Oxygen Concentrator |  |

- Compact & lightweight.
- High efficiency, low noise. Excellent heat dissipation.
- 93% O<sub>2</sub> concentration.
- High adsorption sieve bed.
- Includes AC/DC Car Charger.
- Conforms to all applicable FAA acceptance criteria for POC carriage and use on board aircraft.

## **FEATURES**

| Size                    | 8.7" x 3.4" x 6.3"  |  |
|-------------------------|---|--|
| Weight                  | 4.37lbs   |  |
| Warm up time            | 2 minutes   |  |
| Settings                | Seven Flow Settings (1-7)   |  |
| Oxygen Concentrator     | 87%-96%   |  |
| AC Power Supply:        | 100-240v, 50-60hz (auto sensing to allow for worldwide use)                           |  |
| DC Power Supply:        | Input: 11-16VDC, Output: 19V 6.3a   |  |
| Rechargeable Battery    | 14.4VDC   |  |
| Battery Rated Capacity: | 6.7Ah   |  |
| Display Screen          | 2.8"  |  |
| Battery                 | 6700mAh, duration up to 6 hours on setting 1.<br>Approx. 4 hours charge time with AC. |  |
| Operation               | Simple control functions and easy-to-read LCD display                                 |  |
| Use                     | Designed for 24/7 use   |  |

| Setting | Battery Life |
|---------|--------------|
| 1       | 6 h          |
| 2       | 5.2 h        |
| 3       | 3.6 h        |
| 4       | 2.8 h        |
| 5       | 2.3 h        |
| 6       | 1.9 h        |
| 7       | 1.4 h        |

## **WARRANTY**

Three years, excluding battery and sieve bed

V8: 3/27/24