

AURA-LT PORTABLE AEROSOL MACHINE

W/ 1 REUSABLE NEB KIT
AND 1 DISPOSABLE NEB KIT



ITEM #	DESCRIPTION
PCN2318B	Aura-LT Portable Aerosol Machine w/ One Reusable Neb Kit and One Disposable Neb Kit

- 120V, 60Hz
- 50W ±10%
- Flow 8-10L
- Product weight 3 lb 1 ounce
- 6.6" x 5.6" x 3.7"

THREE YEAR LIMITED WARRANTY

Warranted to be free from defects in material, and workmanship for a period of 3 years from the date of purchase. The warranty covers only the compressor. The replaceable components like nebulizer, mouthpiece, air hose, and filters are not included. Opening or altering the device invalidates the warranty. The warranty does not cover damage caused by improper handling, accidents, or non-compliance with the operating instructions.