



TRANSPORTABLE OXYGEN CONCENTRATOR

Finally, a Transportable Oxygen Concentrator Built for the Realities of Care.

The TOC3 was designed to fill the industry's biggest gap: a true hospital-to-home solution that supports effortless discharge and enables worry-free, 24/7 oxygen use—whether your patient is flying, walking, or sleeping. With continuous flow up to 3 LPM, pulse settings up to 10, and best-in-class battery performance, the TOC3 delivers the freedom and reliability today's oxygen users deserve.



SPECIFICATIONS

Settings	3 LPM Continuous Flow & 10 Pulse Dose Delivery Settings
Size	11.1" x 6.7" x 15"
Product Weight	17.2 lbs (with battery)
Oxygen Concentration	90% - 3% / + 5.5% (at all settings)
Respiratory Rate	10 to 40 BPM
Display Screen	2.8" Color LCD
Operation	Simple Controls, Easy-to-Read LCD Display
Use	Up to 24/7 Use
AC Input	100 to 240V, 50 to 60 Hz
Rechargeable Battery	14.4V DC
Battery Duration	Up to 6.6 Hours of Battery Life, Up to 4.5 Hours Charge Time
Accessories	Bag & Cart Included for Convenient Transport
HCPCS	E1390 and E1392

Continuous Flow	Battery Life
Setting 1.0 L	3.1 h
Setting 1.5 L	2.5 h
Setting 2.0 L	1.9 h
Setting 2.5 L	1.4 h
Setting 3.0 L	1.2 h

Pulse Dose	Battery Life
Setting 1	6.6 h
Setting 2	4.9 h
Setting 3	3.9 h
Setting 4	3.6 h
Setting 5	3.3 h
Setting 6	3.0 h
Setting 7	2.6 h
Setting 8	2.3 h
Setting 9	2.0 h
Setting 10	1.8 h

ACCESSORIES

Description	Item #
Cart	TOCCRT
Carry Bag	TOCCB
Battery	TOCBY-2
AC Power Supply	TOCAC
DC Power Supply	P2DCA-1
Small Intake Filter	TOCIF-1
Large Intake Filter	TOCIF-2

WARRANTY

Covered by our industry leading Rhythm Resolve Warranty program.

3 Years: TOC3

2 Years: Sieve Bed

1 Year: Battery and Accessories



Shown with
Carry Bag
and Cart
(Included)