

RALLY LITE ALUMINUM 3 WHEEL WALKER W/ TOTE

250 lb Weight Capacity



ITEM #	DESCRIPTION
729	Rally Lite Aluminum 3 Wheel Walker with Tote

FEATURES

- Lightweight, 3 wheel rollator with 8" wheels and locking loop brakes
- Telescopic side frame feature makes walker more compact for travel and storage
- Tote bag included

SPECIFICATIONS

Product Weight: 10 lbs
Handle Height Range: 31" - 38"
Product Dimensions: 25" W x 24" L
Weight Capacity: 250 lbs

LIMITED LIFETIME WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service. This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced. For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.