

ARRAY HD HEAVY DUTY EXTRA WIDE K7 WHEELCHAIR

600 lb Weight Capacity



ITEM #	DESCRIPTION
M170-5-226	26" Detachable Desk Length Arms
M170-5-228	28" Detachable Desk Length Arms
M170-5-230	30" Detachable Desk Length Arms
PATLR	Articulating Leg Rests for K1, K2, K3 and K7 Wheelchairs
PELR-HD	ELR Rigging for HD K7 Wheelchairs
PSFR-HD	SF Rigging for HD K7 Wheelchairs

^{*}Front Riggings sold separately.

FEATURES

- Carbon steel frame in silver vein, scratch-proof finish.
- Detachable, reversible padded desk length arms.
- Swing away footrests or elevating leg rests.
- Height adjustable 8" front casters, 3 positions.
- Urethane tires, composite wheels for low maintenance, durability.
- Standard dual axle transitions easily to hemi-level.
- Easy-to-clean, attractive vinyl upholstery.
- Folds to 13" wide for convenient travel and storage.

SPECIFICATIONS	26"	28″	30"
Handle Height	36"	36"	36″
Seat to Floor Height	17.5" - 19.5"	17.5" - 19.5"	17.5" - 19.5"
Seat Dimensions (WxD)	26" x 20"	28" x 20"	30" x 20"
Product Weight	58 lbs	61 lbs	63 lbs
Weight Capacity	600 lbs	600 lbs	600 lbs
Width Open/Closed	32"/13"	34"/13"	36"/13"
HCPC Code	K0007	K0007	K0007

V4: 01/18/22



WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.