

ARRAY K2 WHEELCHAIR

300 lb Weight Capacity



| ITEM # | DESCRIPTION |
|---------------|---|
| M170-5-216SF | 16" Detachable Desk Length Arms & Swing Away Foot Rests |
| M170-5-216ELR | 16" Detachable Desk Length Arms & Elevating Leg Rests |
| M170-5-218SF | 18" Detachable Desk Length Arms & Swing Away Foot Rests |
| M170-5-218ELR | 18" Detachable Desk Length Arms & Elevating Leg Rests |
| M170-5-220SF | 20" Detachable Desk Length Arms & Swing Away Foot Rests |
| M170-5-220ELR | 20" Detachable Desk Length Arms & Elevating Leg Rests |

FEATURES

- Carbon steel frame in silver vein, scratch-proof finish.
- Detachable, reversible padded desk arms.
- Swing away footrests or elevating leg rests.
- Height adjustable 8" front casters, 3 positions.
- Urethane tires, composite wheels for low maintenance, durability.
- Standard dual axle transitions easily to hemi-level.
- Easy-to-clean, attractive vinyl upholstery.
- Folds to 12.5" wide for convenient travel and storage.

| SPECIFICATIONS | 16" | 18" | 20" |
|-----------------------|---------------|---------------|---------------|
| Handle Height | 36" | 36" | 36" |
| Seat to Floor Height | 17.5" - 19.5" | 17.5" - 19.5" | 17.5" - 19.5" |
| Seat Dimensions (WxD) | 16" x 16" | 18" x 16" | 20" x 16" |
| Product Weight | 40.75 lbs | 40.75 lbs | 40.75 lbs |
| Weight Capacity | 300 lbs | 300 lbs | 300 lbs |
| Width Open/Closed | 24"/12.5" | 26"/12.5" | 28"/12.5" |
| HCPC Code | K0001/K0002 | K0001/K0002 | K0001/K0002 |

V6: 1/13/20

ARRAY
K2 WHEELCHAIR
300 lb Weight Capacity



LIMITED LIFETIME WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.