

# COMMANDO K4 WHEELCHAIR

## 300 lb Weight Capacity



ITEM #	DESCRIPTION
M170-3-416SF	16" Adjustable Height Desk Length Arms & Swing Away Foot Rests
M170-3-416ELR	16" Adjustable Height Desk Length Arms & Elevating Leg Rests
M170-3-418SF	18" Adjustable Height Desk Length Arms & Swing Away Foot Rests
M170-3-418ELR	18" Adjustable Height Desk Length Arms & Elevating Leg Rests
M170-3-420SF	20" Adjustable Height Desk Length Arms & Swing Away Foot Rests
M170-3-420ELR	20" Adjustable Height Desk Length Arms & Elevating Leg Rests

### FEATURES

- Carbon steel frame in silver vein, scratch-proof finish.
- Detachable, reversible padded desk arms.
- Swing away footrests or elevating leg rests.
- Height adjustable arms
- Height adjustable 8" front casters, 3 positions.
- Urethane tires, composite wheels for low maintenance, durability.
- Quick release, rear wheels.
- Folds to 12" wide for convenient travel and storage.
- Velcro patch under seat for additional support and comfort.

SPECIFICATIONS	16"	18"	20"
Handle Height	36"	36"	36"
Seat to Floor Height	17.5" - 19.5"	17.5" - 19.5"	17.5" - 19.5"
Seat Dimensions (WxD)	16" x 16"	18" x 16"	20" x 16"
Product Weight	33.5lbs*	33.5 lbs*	33.5 lbs*
Weight Capacity	300 lbs	300 lbs	300 lbs
Width Open/Closed	23"/12"	25"/12"	27"/12"
HCPC Code	K0004	K0004	K0004

\*without front rigging

# COMMANDO K4 WHEELCHAIR

## 300 lb Weight Capacity

### LIMITED LIFETIME WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.