

FLOW K1 WHEELCHAIR

300 lb Weight Capacity



ITEM #	DESCRIPTION
M170-1-116SF	16" Flip Back Desk Length Arms & Swing Away Foot Rests
M170-1-116ELR	16" Flip Back Desk Length Arms & Elevating Leg Rests
M170-1-118SF	18" Flip Back Desk Length Arms & Swing Away Foot Rests
M170-1-118ELR	18" Flip Back Desk Length Arms & Elevating Leg Rests
M170-1-120SF	20" Flip Back Desk Length Arms & Swing Away Foot Rests
M170-1-120ELR	20" Flip Back Desk Length Arms & Elevating Leg Rests

FEATURES

- Carbon steel frame in scratch-proof finish.
- Swing away footrests or elevating leg rests.
- Attractive nylon upholstery is durable and comfortable.
- Urethane tires, composite wheels for low maintenance, durability.
- Folds to 11" wide for convenient travel and storage.

SPECIFICATIONS	16"	18"	20"
Handle Height	36"	36"	36"
Seat to Floor Height	19.5"	19.5"	19.5"
Seat Dimensions (WxD)	16" x 16"	18" x 16"	20" x 16"
Product Weight	39.25 lbs	39.25 lbs	39.25 lbs
Weight Capacity	300 lbs	300 lbs	300 lbs
Width Open/Closed	22"/11"	24"/11"	26"/11"
HCPC Code	K0001	K0001	K0001

V6: 09/07/22

FLOW K1 WHEELCHAIR

300 lb Weight Capacity

LIMITED LIFETIME WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.