

# FULL ELECTRIC BED

## 450 lb Weight Capacity



ITEM #	OPTIONS	HCPC CODE
H9500FULL	Full Electric Bed Only	E0297
H9500FF	Full Electric Bed w/ Full Rails w/o Mattress	E0266
H9500FH	Full Electric Bed w/ Half Rails w/o Mattress	E0266
H9500FFM	Full Electric Bed w/ Full Rails, w/ Mattress	E0265
H9500FPKG	Full Electric Bed w/ Half Rails, w Mattress	E0265

### FEATURES

- UL approved motor.
- Frame: 11ga (.120) hot rolled steel (roll formed).
- Zinc plated steel wire and coils.
- Double powder coated frame.
- Requires a 9V battery (not included).
- Accommodates a 36" x 80" mattress.

### SPECIFICATIONS

SPECIFICATIONS	H9500
Bed Height Range	15" - 24"
Bed Length	87.9"
Bed Width	36.2"
Product Weight	169 lbs
Shipping Weight	182 lbs
Weight Capacity	450 lbs

V3: 1/12/21

# FULL ELECTRIC BED

## 450 lb Weight Capacity



Raise Head — Lower Head  
Raise Foot — Lower Foot  
Raise Head & Foot Together — Lower Head & Foot Together

FULL PENDANT CONTROL

### WARRANTY

Limited Lifetime warranty against manufacturing defects on welds.  
5 year warranty on Frame and Motor. 1 year warranty on all other parts and components.

This warranty does not cover product failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. It is recommended but not required that the product be returned to the Rhythm Healthcare, LLC dealer through whom it was originally purchased. However, all Rhythm Healthcare, LLC Dealers are qualified to assist you in obtaining warranty service. If the product is to be returned to Rhythm Healthcare, LLC, prior authorization will be required. In the event of a defect covered by this warranty, Rhythm Healthcare, LLC will determine if the product will be repaired or replaced.

This Rhythm Healthcare, LLC product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have questions about your Rhythm Healthcare, LLC product or this warranty, please email us at [contactus@rhythmhc](mailto:contactus@rhythmhc).