

# DELUXE RECLINING DESK LENGTH ARMCHAIR

**300 lb Weight Capacity**



ITEM #	DESCRIPTION
L9016	16" Deluxe Reclining Desk Length Armchair
L9018	18" Deluxe Reclining Desk Length Armchair
L9020	20" Deluxe Reclining Desk Length Armchair

## FEATURES

- Removable desk length arms
- Comes standard with ELRS
- Vinyl upholstery
- Quick release rear wheels
- Available Color: Hammertone Grey

SPECIFICATIONS	16"	18"	20"
Handle Height	36"	36"	36"
Seat to Floor Height	20"	20"	20"
Seat Dimensions (WxD)	16" x 16"	18" x 16"	20" x 16"
Product Weight	53 lbs	54.5 lbs	55.6 lbs
Weight Capacity	300 lbs	300 lbs	300 lbs
Product Dimensions	26" x 36"	26" x 36"	26" x 36"

## WARRANTY

Rhythm Healthcare, LLC offers a limited 1 year warranty from defects in materials and craftsmanship of the product to the original purchaser. This warranty does not cover normal wear and tear or misuse of the product. Rhythm Healthcare, LLC reserves the right to determine if a warranted product should be repaired or replaced.

V2: 9/5/24