

ROYAL SERIES ROLLATORS



807BB
(Blueberry)

ITEM	COLORS
807S - Royal Steel	807SBB (Blueberry) 807SGT (Garnet)
807 Royal Aluminum	807AS (Amber Sunrise) 807BL (Blue) 807BB (Blueberry) 807CR (Candy Apple Red) 807RD (Red) 807SG (Slate Grey)
807U Royal Universal	807ULB (Laser Blue) 807ULR (Laser Red) 807ULP (Laser Purple)
812U Royal Deluxe Universal	812UAQ (Aqua) 812ULB (Laser Blue) 812ULG (Laser Green) 812ULR (Laser Red)
850U Royal Bariatric Universal	850ULB (Laser Blue) 850ULR (Laser Red)

NEW COLORS!



807BB
(Blueberry)



807SBB
(Blueberry)



807SGT
(Garnet)



807U
(Blueberry)



807CR
(Candy Apple Red)



807SG
(Slate Grey)



807BB
(Blueberry)



807CR
(Candy Apple Red)



807SG
(Slate Grey)



807CR
(Candy Apple Red)



807SG
(Slate Grey)

ROYAL SERIES ROLLATORS



SPECIFICATIONS	807	807S
Weight	14 lbs	17 lbs
Handle Height	30" - 37"	32" - 38"
Height Adjustment	Handles only	Handles only
Seat to Floor Height	20.25"	22"
Seat Dimensions	17.5" W x 13" D	17.5" W x 13" D
Unit Dimensions	23.25" W x 26" L	23.25" W x 26" L
Frame	Aluminum Round	Steel Round
Wheels	6"	6"
Weight Capacity	300 lbs	300 lbs

SPECIFICATIONS	807U	812U	850U
Weight	15 lbs	17 lbs	20 lbs
Handle Height	30.5" - 38"	1.5" - 39"	31.5" - 40"
Height Adjustment	Handle & Seat	Handle & Seat	Handle & Seat
Seat to Floor Height	18.5" - 21.5"	19.5" - 22.5"	18.5" - 22"
Seat Dimensions	18.5" W X 14" D	14.5" W X 14.5" D	18" W X 15" D
Unit Dimensions	23.25" W X 26" L	24.5" W X 27.5" L	27" W x 28.5" L
Frame	Aluminum Oval	Aluminum Oval	Aluminum Oval
Wheels	6"	8"	8"
Weight Capacity	300 lbs	300 lbs	450 lbs

WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this **Limited Lifetime Warranty** to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.