

# TWO BUTTON FOLDING WALKER

## 300 lb Weight Capacity



**Adjusts from junior to adult:**  
28.25" to 38.25 (w/out wheels)  
28.25" to 38.5 (w/ wheels)

### NEW!

Walker w/ 5" Wheels  
is now available  
assembled in a 1 pack.  
**USE ITEM # P1307.**

ITEM#	DESCRIPTION
P1300	Two Button Folding Walker (Without Wheels)
P1305 (Assembly Required, 4/case)	Two Button Folding Walker (With 5" Wheels)
P1307 (Assembled, 1/case)	Two Button Folding Walker (With 5" Wheels)

## FEATURES

- Universal Height Adjustment from Junior to Tall Adult

SPECIFICATIONS	P1300	P1305	P1307 (Assembled)
Weight	6 lbs	6.5 lbs	6.5 lbs
Height	28.25" - 38.25"	28.5" - 38.5"	28.5" - 38.5"
Unit Dimensions (WxD)	22.25" X 17.25"	24.5" X 19"	24.5" X 19"
Wheels	x	✓	✓
Assembled	x	x	✓

## WARRANTY

Rhythm Healthcare, LLC offers a limited 1 year warranty from defects in materials and craftsmanship of the product to the original purchaser. This warranty does not cover normal wear and tear or misuse of the product. Rhythm Healthcare, LLC reserves the right to determine if a warranted product should be repaired or replaced

V7: 9/12/24