

# CADENCE K3 WHEELCHAIR

300 lb Weight Capacity



ITEM #	DESCRIPTION
M170-1-316SF	16" Detachable Desk Length Flip Back Arms & Swing Away Foot Rests
M170-1-316ELR	16" Detachable Desk Length Flip Back Arms & Elevating Leg Rests
M170-1-318SF	18" Detachable Desk Length Flip Back Arms & Swing Away Foot Rests
M170-1-318ELR	18" Detachable Desk Length Flip Back Arms & Elevating Leg Rests
M170-1-320SF	20" Detachable Desk Length Flip Back Arms & Swing Away Foot Rests
M170-1-320ELR	20" Detachable Desk Length Flip Back Arms & Elevating Leg Rests
M170-3-316SF	16" Adjustable Height Desk Arms & Swing Away Foot Rests
M170-3-316ELR	16" Adjustable Height Desk Arms & Elevating Leg Rests
M170-3-318SF	18" Adjustable Height Desk Arms & Swing Away Foot Rests
M170-3-318ELR	18" Adjustable Height Desk Arms & Elevating Leg Rests
M170-3-320SF	20" Adjustable Height Desk Arms & Swing Away Foot Rests
M170-3-320ELR	20" Adjustable Height Desk Arms & Elevating Leg Rests

## FEATURES

- Aluminum frame with steel cross brace.
- Flip back detachable, padded desk arms for easy self transfer.
- Swing away footrests or elevating leg rests.
- Height adjustable 8" front casters, 3 positions.
- Urethane tires, composite wheels for low maintenance, durability.
- Standard dual axle transitions easily to hemi-level.
- Folds to 12" wide for convenient travel and storage.
- Attractive nylon upholstery is durable and comfortable.
- Optional add on: height adjustable arms.

HEIGHT  
ADJUSTABLE  
ARMS



V10: 3/03/21

# CADENCE K3 WHEELCHAIR

300 lb Weight Capacity



SPECIFICATIONS	16"	18"	20"
Handle Height	36"	36"	36"
Seat to Floor Height	17.5" - 19.5"	17.5" - 19.5"	17.5" - 19.5"
Seat Dimensions (WxD)	16" x 16"	18" x 16"	20" x 16"
Product Weight	35 lbs*	35 lbs*	35 lbs*
Weight Capacity	300 lbs	300 lbs	300 lbs
Width Open/Closed	23"/12"	25"/12"	27"/12"
HCPC Code	K0003	K0003	K0003

\*Without front rigging

## LIMITED LIFETIME WARRANTY

Lifestyle Mobility Aids warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC - Lifestyle Mobility Aids dealer.