



RHYTHM
HEALTHCARE

OXYGEN REGULATORS

Our Oxygen Regulators feature a 1.5" easy-to-read gauge and an internal pressure relief valve for safety. With our mini regulators' compact size, they easily fit in portable oxygen bags. Mini regulators provide equivalent performance to our standard size regulators.

- Aluminum body with brass core
- 1.5" pressure gauge
- Easy-to-read flow rate window
- Easy-to-operate single selector knob
- Lightweight and durable design
- Complies to CGA standards

MINI OXYGEN REGULATORS

Mini Oxygen Regulator, 0 - 8 LPM , Barb Outlet, CGA-870



ITEM#: OR108B

Specifications

- Barb Outlet
- Mini Standard Flow
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved
- Connection Type: CGA-870
- Flow Control Type: Click Style
- Flow Rate: 0 - 8 LPM

Mini Oxygen Regulator, 0 - 15 LPM , Barb Outlet, CGA-870



ITEM#: OR1015B

Specifications

- Barb Outlet
- Mini Standard Flow
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved
- Connection Type: CGA-870
- Flow Control Type: Click Style
- Flow Rate: 0 - 15 LPM

STANDARD OXYGEN REGULATORS

Oxygen Regulator, 0 - 8 LPM , Barb Outlet, CGA-870

ITEM#: OR208B



Specifications

- Barb Outlet
- Standard Flow
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved
- Connection Type: CGA-870
- Flow Control Type: Click Style
- Flow Rate: 0 - 8 LPM

Oxygen Regulator, 0 - 15 LPM , Barb Outlet, CGA-870

ITEM#: OR2015B



Specifications

- Barb Outlet
- Standard Flow
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved
- Connection Type: CGA-870
- Flow Control Type: Click Style
- Flow Rate: 0 - 15 LPM

Oxygen Regulator, 0 - 8 LPM , Diss Outlet, CGA-540

ITEM#: OR608D



Specifications

- Diss Outlet
- Standard Flow
- Inlet Pressure: 300 - 3,000 PSI
- Outlet pressure: 50 PSI
- UL Approved
- Connection Type: CGA-540
- Flow Control Type: Click Style
- Flow Rate: 0 - 8 LPM