

# ARPEGGIO ROLLATOR

## 300 lb Weight Capacity



ITEM #	DESCRIPTION
955CH	Champagne

### FEATURES

- Backrest is height and depth adjustable for optimal comfort.
- Cross-brace frame for added strength and safety.
- Quick and easy to fold up.
- Large removable tote bag for convenient storage included.
- Top of the line leather backrest and seat for comfort and durability.
- 8" Chrome Wheels for indoor/outdoor use.
- Clamp lock holds rollator closed when it's folded.

### SPECIFICATIONS

Weight	17 lbs
Handle Height	32" - 37"
Height Adjustment	Handles Only
Seat to Floor Height	19.75"
Adjustable Backrest	4" Height Adjustment and 1.5" Depth Adjustment
Seat Dimensions	18" x 8.75"
Unit Dimensions	23" W x 28.5" L
Frame	Aluminum
Weight Capacity	300 lbs

### WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this Limited Lifetime Warranty to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.