

19" ALUMINUM COMPANION CHAIR

300 lb Weight Capacity



ITEM #	DESCRIPTION
L2019BL	19" Aluminum Companion Chair (Blue)
L2019BG	19" Aluminum Companion Chair (Burgundy)

FEATURES

- Adjustable swing away footrests
- 8" Wheels
- Nylon Upholstery
- Convenient fold down back for storage and transport
- Standard rear wheel locks

SPECIFICATIONS

Weight	20 lbs
Handle Height	40"
Seat to Floor Height	19"
Seat Dimensions	19"W x 16"D
Unit Dimensions	22.5" x 24"
Weight Capacity	300 lbs

LIMITED LIFETIME WARRANTY

Limited Lifetime Warranty against manufacturing defects. This warranty does not cover product failure due to misuse, negligence, improper storage or handling, improper operation or unauthorized modifications. Normal wear and tear on non-durable components, such as rubber accessories and casters are not covered under this warranty. These items are subject to normal wear and need periodic replacement. For Warranty Service, it is recommended but not required, that the product be returned to the Rhythm Healthcare dealer through whom it was originally purchased. However, all Rhythm Healthcare dealers are qualified to assist you in obtaining warranty service. If the product is to be returned to Rhythm Healthcare prior authorization will be required. In the event of a defect covered by this warranty, Rhythm Healthcare will determine if the product will be repaired or replaced.