

# HEIGHT ADJUSTABLE ROLLATOR - 6" CASTORS

## 300 lb Weight Capacity



ITEM #	DESCRIPTION
807ADJBK	Height Adjustable Rollator (Black)

### FEATURES

- Universal seat height adjustment from 18" to 22" (1" increments)
- Comfortable, padded seat
- Padded backrest, which is removable and hinged, can be folded up or down.
- 6" casters ideal for indoor/outdoor use
- Height adjustable & easy to grip handles

### SPECIFICATIONS

Product Weight:	14 lbs
Frame:	Aluminum
Seat Dimensions:	13" (W) x 12" (D)
Handle Height:	29.5" - 38"
Seat to Floor Height:	18" - 22"
Casters:	6"
Width between Arms:	17"
Weight Capacity:	300 lbs

### WARRANTY

Rhythm Healthcare, LLC offers a limited 1 year warranty from defects in materials and craftsmanship of the product to the original purchaser. This warranty does not cover normal wear and tear or misuse of the product. Rhythm Healthcare, LLC reserves the right to determine if a warranted product should be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.

V2: 8/4/24