

REGAL - BARIATRIC ALUMINUM 4 WHEEL ROLLATOR WITH UNIVERSAL HEIGHT ADJUSTMENT



ITEM	COLORS
850U Royal Bariatric Universal	850ULB (Laser Blue) 850ULR (Laser Red)

SPECIFICATIONS 850U

Weight	20 lbs
Handle Height	31.5" - 40"
Between Handles	21.5"
Height Adjustment	Handle & Seat
Seat to Floor Height	18.5" - 22"
Seat Dimensions	18" x 15"
Unit Dimensions	27" W x 28.5" L
Frame	Aluminum Oval
Wheels	8"
Weight Capacity	450 lbs

WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this **Limited Lifetime Warranty** to affirm our confidence that the materials and workmanship of our products will result in years of dependable service. This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.