

PEDIATRIC COMPRESSOR NEBULIZER



ITEM #S:
CCN2318BL,
CCN2318PNK



SPECIFICATIONS	
Operating Pressure	13 PSI
Particle Size	0.5 μ m to 5 μ m
Weight	3.1 lbs
Electrical Requirements	120 VAC 60Hz
Liter Flow	8 lpm
Method of Nebulization	Piston Compressor
Operating Pressure	13 PSI
HCPCS Code	E0570-1
Includes	Disposable nebulizer kit, pediatric aerosol mask, and replacement filter
Available Colors	Blue (CCN2318BL), Pink (CCN2318PNK)
Warranty	5 Years



Federal law restricts this device to sale by, or on the order of, a licensed healthcare practitioner



Kid-Friendly: Makes treatments more engaging and encourages consistent use.



Efficient & Quiet Pump: Powerful piston pump delivers dependable performance while maintaining quiet operation.



Complete Pediatric Kit: Includes a disposable nebulizer kit, pediatric aerosol mask, and 5 replacement filters.