

INTRODUCING THE 50 PSI COMPRESSOR



The Rhythm Healthcare 50psi Compressor is designed for patients who require both humidification and medication delivery through nebulization. With its powerful compressor, it provides up to 50 psi of maximum pressure, ensuring fast and efficient medication delivery. Engineered for quiet operation, the Rhythm Healthcare 50 psi Compressor is portable, easy to clean, and requires minimal maintenance, making it a reliable choice for patients and caregivers alike.



- » Provides efficient, continuous high-pressure performance
- » Ultra-powerful compressor for both humidification and high flow, high-pressure nebulization of medications, including viscous formulations
- » Easy to Clean and maintain
- » Delivers up to 50psi at maximum pressure

ORDERING INFORMATION

ITEM #	Description
CMP501	50 psi Compressor
CMP501-2	Intake Filter
CMP501-3	Filter Cap

TECHNICAL SPECIFICATIONS

Size	12.2 x 7.7 x 15.3 in
Weight	17.2 lbs
Operating Flow Range	11-34 LPM
Operating Pressure Range	7-50 psi
Operating Temperature Range	50°F - 104°F

