

SOFTCAIRE EC (WITH FOAM BASE)

Weight Capacity: 350 lbs



ITEM #	DESCRIPTION
H80513	SoftCaire EC (With Foam Base)

Quick Connector with Cap

- Ensures Mattress remains inflated during power outage or patient transport.

CPR function

- Allows rapid deflation in case of medical emergencies.
- Designed for one handed operation.
- Positioned at head end.

Pillow Function

- Head zone comprised of three fully inflated static air cells to provide patient comfort.

Foam Base

System includes 17 individual air cells over a 3" foam base as an effective pressure support base.

SPECIFICATIONS

Pump:

Weight Setting Control:	Dial
Dimension:	10.83"L x 3.94"W x 6.10" H
Weight:	5.5 lbs
Air Output:	8 liter/min
Case Material:	Flame Retardant Abs
Pressure Range:	30-60 mmHg
Cycle Time:	10 min
Power Supply:	120V/60 Hz
Power Cord:	14' 9"

Mattress:

Inflated Dimension:	80" L x 35.43" W x 4.92" H
Alternating:	A/B Alternating
Structure:	17 x 5" Air Bladders
Cover Material:	Nylon/PU Quilted with Zipper
Cell Material:	Nylon/PVC
Base Material:	Nylon/PVC
Weight Capacity:	350 lbs
Flame Retardant Standards:	CA117/CFR 1632

SOFTCAIRE EC
(WITH FOAM BASE)
Weight Capacity: 350 lbs



18 MONTH LIMITED WARRANTY

Your Rhythm product is warranted to be free of defects in materials and workmanship for 18 Months to the original consumer purchaser. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, which are subject to normal wear and need periodic replacement.

If the product is to be returned to Rhythm Healthcare, prior authorization will be required. This product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have a question about your Rhythm Healthcare product or this warranty, please email us at contactus@rhythmhc.com