

SOFTCAIRE BASIC

Weight Capacity: 350 lbs



ITEM #	DESCRIPTION
H80521	SoftCaire Basic

The SoftCaire Basic Group 2 Support Surface is a combination therapy Alternating Pressure and Low Air Loss mattress system for cost-effective, efficient and effective pressure redistribution treatment.

Cover

- Top Cover is Low Shear, Water-Resistant, Vapor-Permeable, and Anti-Microbial

Quick Connector with Cap

- Ensures Mattress Remains Inflated during power outage or patient transport

CPR Function

- Allows for One-Handed, Rapid Deflation in case of medical emergencies or to facilitate rapid disassembly

Pillow Function

- Head Zone Comprised of Three Fully Inflated Static Air Cells to provide patient comfort

Heel Function

- Provides Single Cell Deflation to Float Vulnerable Heels

Spring-Loaded, Anti-Kink Tubes

- Enhances experience and reduces service calls

SPECIFICATIONS

Pump

Weight Setting Control:	Dial
Dimension:	10.83"L x 3.94"W x 6.10" H
Weight:	5.5 lbs
Air Output:	8 liter/min
Case Material:	Flame Retardant Abs
Pressure Range:	30-60 mmHg
Cycle Time:	10 min
Power Supply:	120V/60 Hz
Power Cord Length:	14' 9"

Mattress

Inflated Dimension:	80"L x 36"W x 8"H
Alternating:	A/B Alternating
Structure:	20 x 8" Air Bladders
Air cells:	Nylon/ PVC with low air loss
Base:	70D Nylon/PVC base cover with 3 anchoring straps on each side
Maximum Weight Capacity:	350 lbs
Flame Retardant Standards:	CA117

SOFTCAIRE BASIC

Weight Capacity: 350 lbs



18 MONTH LIMITED WARRANTY

Your Rhythm product is warranted to be free of defects in materials and workmanship for eighteen (18) months to the original consumer purchaser. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, which are subject to normal wear and need periodic replacement.

If the product is to be returned to Rhythm Healthcare, prior authorization will be required. This product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have a question about your Rhythm Healthcare product or this warranty, please email us at contactus@rhythmhc.com