

THE RHYTHM PULSEPROTECT MATTRESS



ITEM #	DESCRIPTION
H81012	The Rhythm PulseProtect Mattress

- Provides superior pressure redistribution and is ideal for low to moderate risk patients.
- Channel Cut High Density Polyurethane Foam. Individual cells offer superior pressure redistribution.
- Nylon top cover is vapor permeable which reduces perspiration and provides a low shear/friction surface that is anti-microbial/anti-bacterial and fluid resistant.
- 360° Zipper for easy removal, cleaning or replacement.
- Anti-Contamination Flap covers the zipper reducing the possibility of fluids damaging the foam core.
- Bottom Cover: Triple laminate vinyl; fluid resistant; anti-microbial/anti-bacterial; fire retardant.
- Fire sleeve meets Federal Fire Code 16 CFR part 1633.
- Weight Capacity: 350 lbs.
- 36" x 80" x 6"

WARRANTY

Limited Warranty: Foam 2 Years, Cover 1 Year

V2: 2/11/25