

SEMI-ELECTRIC LOW HOMECARE BED



ITEM #	DESCRIPTION
H9200SF	Semi-Electric Low Homecare Bed with Full Side Rails
H9200SH	Semi-Electric Low Homecare Bed with Half Side Rails

SPECIFICATIONS

Overall dimension	35.63" W x 87.2"L
Sleep surface	34.65" W x 79.73" L
Mattress:	80" mattress required
Lowest height (inner castor setting)	15.35" H
Lowest height (outer castor setting)	23.25" H
Motor:	1x
Castor	2 x 3" castor w/locking and 2 x 3" castor w/o locking
Hand Crank & Gear Box:	Plastic
Weight Capacity	450 lbs

WARRANTY

Original Purchaser. Limited Lifetime on welds, 5 year on frame, motor and pendant. 1 Year on all other parts and components. During the warranty period, Rhythm Healthcare will determine if a defective item will be repaired or replaced.