

BARIATRIC EXTRA WIDE TWO BUTTON FOLDING WALKER

500 lb Weight Capacity



Adjusts from:
32" to 39"

ITEM#	DESCRIPTION
P1500	Without Wheels
P1505	With 5" Wheels

FEATURES

- For bariatric use
- Wider and deeper frame design

SPECIFICATIONS	P1500	P1505
Weight	7.5 lbs	9 lbs
Height adjustable	32" to 39"	32" to 39"
Unit Dimensions (WxD)	24.5" W x 20" L	25.75" W x 21.5" L

WARRANTY

Rhythm Healthcare, LLC offers a limited 1 year warranty from defects in materials and craftsmanship of the product to the original purchaser. This warranty does not cover normal wear and tear or misuse of the product. Rhythm Healthcare, LLC reserves the right to determine if a warranted product should be repaired or replaced

V2: 9/15/23