

INCREASED MOBILITY FOR ELEVATED OXYGEN NEEDS



ITEM #: P2-E6



P2-E6 PORTABLE OXYGEN CONCENTRATOR

More oxygen support in a lightweight design.

Meet higher clinical demands without compromising patient mobility. The P2-E6 gives providers a flexible portable solution for oxygen users with increased oxygen needs. Built for everyday mobility and simple operation, it offers a practical step up from standard five-setting POCs without adding unnecessary size or complexity.



6 Pulse Settings

- ✓ **Expanded Clinical Capacity:** The P2-E6 delivers up to 1,200 ml/min across 6 flow settings, allowing providers to comfortably accommodate users with elevated oxygen needs.
- 🍃 **Lightweight, Extended Range:** Weighing just 4.37 lbs with up to 5.4 hours of battery life, it delivers the ideal balance of lightweight portability and long-lasting freedom.
- 🛡️ **Unparalleled Support:** Covered by the Rhythm Resolve™ Warranty, guaranteeing fast replacement instead of repair throughout the industry-leading 5-year warranty period.

P2-E6 PRODUCT INFORMATION



ITEM #: P2-E6

SPECIFICATIONS

Settings	Six Flow Settings (1-6)
Dimensions	8.7" L x 3.4" W x 6.3" H
Product Weight	4.37 lbs
Warm up time	2 Minutes
Oxygen Concentration	87% - 96%
Respiratory Rate	10 to 40 breaths/min
Maximum Oxygen Output	1,200 ml/min
Sound Level	≤36 dBA (on Setting 2)
Auto-Pulse Technology	Auto-Pulse after 60 Seconds
Auto-Adjustable Trigger Sensitivity	-0.07 to -0.12 cm H ₂ O (5-level sensitivity adjustment)
Alarms	Power Failure, Temperature, Low Oxygen, Battery, Pressure, Sieve Bed Replacement
AC Power Supply	100-240v, 50/60hz; worldwide auto-sensing
DC Power Supply	DC Input: 11 to 16 V DC Output: 19 V 6.3 A
Battery	6700mAh, 14.4V Approx. 4h Charge time with AC
Display Screen	2.8" LCD
Operation	Simple Controls, Easy-to-Read LCD Display
Use	Designed for 24/7 Use
Included Accessories	AC/DC Car Charger, Carry Bag
FAA Compliance	Meets FAA acceptance criteria
HCPCS Code	E1390, E1392
Warranty	Device: 5 Years Sieve Bed: 2 Years Battery: 1 Year

SETTINGS

Breath Rate	1	2	3	4	5	6
	Pulse Volumes (ml)					
10	21	42	63	84	100	120
15	14	28	42	56	66.7	80
20	10.5	21	31.5	42	50	60
25	8.4	16.8	25.2	33.6	40	48
30	7	14	21	28	33.3	40
35	6	12	18	24	28.6	34
40	5.3	10.5	15.8	21	25	30

SETTING	BATTERY LIFE
1	5.4 hr
2	4.8 hr
3	3.3 hr
4	2.5 hr
5	2.1 hr
6	1.7 hr

Rhythm Resolve™ Replacement Warranty



The P2-E6 is covered by our industry-leading Rhythm Resolve™ Replacement Warranty.

For eligible units, the 5-year warranty provides a brand-new replacement rather than a repair in the event of a manufacturer's defect.



SCAN THE QR CODE:

View P2-E6 features, specifications, and resources