



AURA -LT **AEROSOL** DELIVERY SYSTEM

The Aura-LT Aerosol Delivery System is a jet nebulizer that provides fast treatment time at an outstanding value.

The Aura-LT is available in 2 configurations for your convenience. One configuration contains a disposable neb kit and the other contains both a disposable and reuseable neb kit, providing you a choice when serving your patients.



AVAILABLE IN 2 CONFIGURATIONS

PRODUCT SPECIFICATIONS

ITEM #	Description
PCN2318A	Aura-LT Aerosol Machine with Disposable Nebulizer Kit
PCN2318B	Aura-LT Aerosol Machine with One Reusable Nebulizer Kit and One Disposable Nebulizer Kit

TECHNICAL SPECIFICATIONS

Power	120 V AC, 60 Hz
Compressor Flow	> 8 LPM
Weight	3 lbs, 1 oz.
Dimensions	6.6" x 5.6" x 3.7"
Compressor Warranty	5 Years

