

# AURA-NEB PORTABLE NEBULIZER



ITEM #	DESCRIPTION
PCN843	Aura-Neb Portable Nebulizer

- Max Flow 12L/min
- Max pressure 240 kPa
- Nebulization Rate 0.4 ml/min

## THREE YEAR LIMITED WARRANTY

Warranted to be free from defects in material, and workmanship for a period of 3 years from the date of purchase. The warranty covers only the compressor. The replaceable components like nebulizer, mouthpiece, air hose, and filters are not included. Opening or altering the device invalidates the warranty. The warranty does not cover damage caused by improper handling, accidents, or non-compliance with the operating instructions.