

*Wherever you Go
Oxygen Follows*

SPECIFICATIONS:

Model:	P2
Size:	3.3" W x 8.7" D x 6.3" H
Weight:	4.37 lbs
Warm up time:	2 Minutes
Settings:	Five Flow Settings (1-5)
Screen Size:	2.8"
Battery:	6800mAh, duration up to 5 hours on setting 1
Operation:	Simple control functions and easy-to-read LCD display
Use:	
Warranty	Designed for 24/7 use Three years, excluding battery and sieve bed

Provider info:



www.rhythmhc.com
(877) 843-6464
contactus@rhythmhc.com



**P2 PORTABLE
OXYGEN
CONCENTRATOR**





**COMPACT
&
LIGHTWEIGHT**

CONVENIENT OXYGEN 24/7

P2 Portable Oxygen Concentrator eliminates the need for bulky and inconvenient tanks. The P2 is lightweight at only 4.37 lbs.

Features:

- Small & Convenient
- Simple control functions and easy-to-read 2.8" LCD display
- High Absorption Sieve Bed

Includes:

- Nasal Cannula
- Carry Bag with Shoulder Strap
- AC power supply
- DC power supply



Duration up to 5 Hours

Equivalent up to 5 Liters

Lightweight - 4.37 lbs

ACCESSORIES



ITEM #: P2BP-1
Backpack for POC



ITEM #: PAK3
Accessory Kit includes External
Charger, Battery & AC Adapter



ITEM #: P2PBC090
External Battery Charger
(Conveniently charge a backup
battery while the device is in
use.)



ITEM #: P2BY001-2
Battery for POC



ITEM #: P2DCA-1
DC Power Supply



ITEM #: P2ACA-1
AC Power Supply

Purchase accessories at rhythmhc.com



Scan for our catalog.