


1	LM209122	Oxygen outlet fitting
2	LM209122	Flow Meter
3	LM209122	Control panel label
4	LM209122	Hour Meter
5	LM209122	Power switch
6	LM209122	Circuit breaker



7	LM310012	Access Window for the LM5A
8	LM209334	Unit Cabinet for the LM5A
9	LM209065	Caster for the LM5A
10	LM209058	Humidifier for the LM5A
11	LM209136	ST4. 8 x 16 for the LM5A



12	LM209078	Power Cord for the LM5A
13	LM310014	Air Guide Cover for the LM5A