

ROYAL DELUXE UNIVERSAL - ALUMINUM 4 WHEEL ROLLATOR WITH UNIVERSAL HEIGHT ADJUSTMENT



ITEM	COLORS
812U Royal Deluxe Universal	812UAQ (Aqua) 812ULB (Laser Blue) 812ULG (Laser Green) 812ULR (Laser Red)

SPECIFICATIONS

812U

Weight	17 lbs
Handle Height	1.5" - 39"
Height Adjustment	Handle & Seat
Seat to Floor Height	19.5" - 22.5"
Seat Dimensions	14.5" W X 14.5" D
Unit Dimensions	24.5" W X 27.5" L
Frame	Aluminum Oval
Wheels	8"
Weight Capacity	300 lbs

WARRANTY

Rhythm Healthcare, LLC warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer.

We offer this **Limited Lifetime Warranty** to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.